Instructions for use


Talux®

TLXoxyz - Talux without Adapter

TLPoxyz - Talux with Pyramid Adapter

Life Without Limitations®

	3
Ⓔ	Instructions for use.....	6
Ⓕ	Gebrauchsanweisung.....	8
Ⓖ	Guide de fabrication.....	11
Ⓔ	Instrucciones para el uso.....	13
Ⓘ	Istruzioni per l'uso.....	15
Ⓓ	Brugsanvisning.....	17
Ⓔ	Bruksanvisning.....	19
Ⓖ	Gebruiksaanwijzing.....	21
Ⓕ	Instruções de Utilização.....	23

EN - Caution: Össur products and components are designed and tested to ISO 10328. Compatibility and compliance with this standard is achieved only when Össur products and components are used with other recommended Össur or authorized components. If un-usual movement or product wear is detected in a structural part of a prosthesis at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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







ENGLISH

TALUX® SYSTEM (Figure A)

- Indicated for low to moderate impact levels only.
- Build height with Male Pyramid: 176mm (6 15/16").
- Maximum weight 147kg (325 lbs) low activity or 130kg (285 lbs) moderate activity.
- 10 mm (3/8") heel height.
- Available in sizes 23-30.

BENCH ALIGNMENT (Figure B)

- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angels flexion/extension and abduction/adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram.
- The mark denoting the junction of the middle and posterior thirds is located inside the foot cover opening.

DYNAMIC ALIGNMENT

The tapered heel part, compressed tarsal core and the opened dual J-shaped pylons store energy at heel strike, then slowly releasing it to assist tibial progression to midstance. The forward momentum generated by this action results in the toe being loaded for optimum energy release to toe off. The heel to toe action can be influenced by:

- A-P positioning of foot
- Dorsi-Plantarflexion
- Heel Stiffness

HEEL STIFFNESS

Heel too soft

Symptoms

- Foot will come to flat position too early (amputees will feel like sinking into a hole)
- Extra energy is required to climb up over the toe
- Toe will feel too stiff
- Knee may hyperextend

Action

- Add stiffer heel wedge
- Shift socket anterior (or foot posterior)

Heel too hard

Symptoms

- Rapid heel to toe movement
- At heel strike the amputee has poor control of his prosthesis
- Minimal energy return feeling
- Knee may become unstable

Action

- Replace (or remove) wedge with a softer one
- Shift socket posterior (or foot anterior)

HEEL WEDGES, TARSAL CORE AND ACHILLES STRAP (Figure C)

Heel Wedges

- Wedges can be inserted to increase heel stiffness.
- To insert: Slide the wedge between the Talux Tarsal Core and the Achilles Strap retainer until it protrudes evenly on both sides.
- To remove: Grab the exposed end of the wedge with pliers and pull to remove.
- To attach permanently: Trim the wedge to the width of the keel.

- Apply a drop of Loctite 495 to the surface where the wedge meets the foot plate.
- Do not expose Tarsal Core to heat, abrasive materials or excessive moisture as it may cause damage.
- Below are the corresponding stiffness to color:

White	Soft
Red	Medium
Blue	Hard

FLEX-FOOT SOCK

To protect the cover and avoid noise a Flex-Foot sock is fitted over the foot module before installing the cover. The sock should be loose in the heel area.

SHOE HORN (Figure D+E)

When removing and installing the cover use shoehorn to avoid damage to foot or cover.

CAUTION: Don't use the alignment wedges with the male pyramid adapter on the Talux system.

Achilles strap adjustment

To correctly set the achilles strap tension, follow these instructions:

1. Apply Loctite to bolts and fit adapter, finger tight.
2. Loosen up 1/4 turn to free up the hanger
3. Apply tension on strap by pushing the hanger towards top of foot module (**Figure F**)
4. Secure bolts with a wrench to keep hanger in place
5. Check tension on strap.
 - a. No slack should be on the strap when foot is neutral (**Figure G**), check tension by trying to run a piece of paper between strap and tensioner plate, resistance should be felt.
 - b. A too tight strap results in the uppermost carbon plate to lift from the second plate (**Figure H**)
6. Torque to specifications, starting with the upper bolt.
7. Place plastic cap on bolts.

In the case of an achilles strap replacement, the following needs to be carried out:

1. Undo attachment from foot module
2. Using a sharp knife, remove sole material around foot plate bolts to make room for socket
3. Undo foot plate screws with an 8mm socket wrench
4. Clean surfaces if needed with a mild abrasive and/or acetone
5. Put Loctite on foot plate screws and torque the new strap in place with a torque wrench (5 Nm)
6. Follow the instructions above to correctly set strap tension

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Össur Iberia S.A.