ENGLISH

CAUTION

Patients should be instructed to inspect their knee and this product before and after each use.

PHOTO CAPTIONS

- A. Softgood
- B. Posterior Stay
- C. Tibial Cuff and Straps
- D. Middle Adjustable Straps (2)
- E. Femoral Cuff and Straps
- F. Contact Closure Tabs
- G. Locking Buttons and Sizing Indicators

FITTING INSTRUCTIONS

1. Undo all contact closure straps and remove softgood from posterior stay assembly. Apply softgood to the patient, centering softgood knee opening directly over the patella. Ensure pant leg is pulled down taut.
2. Trim softgood width, if necessary, for a custom fit, moving the four contact closure tabs away from the softgood edge as needed.
3. Determine length of posterior stay assembly required for the patient. Product is pre-set to 20" length.
4. Adjust femoral cuff by pressing the blue locking button as shown. Repeat adjustment for the tibial cuff. For best results, place the tibial cuff at or above the calf muscle.
5. Apply posterior stay assembly to the leg. Center the posterior stay along the midline of the leg. Ensure that the softgood extends past the top and bottom ends of the posterior stay assembly. If necessary, trim softgood to length, moving the four contact closure tabs away from softgood edge as needed.
6. Comfortably secure all four straps in order from bottom to top. Ensure that middle adjustable straps are not directly over patella area. If necessary, slide straps up or down for maximum patient comfort. Straps should be positioned at least 2 inches away from the patella.
7. If straps are too long, place an extra contact closure tab at the base of each strap. Attach the strap and trim as desired.

WASHING INSTRUCTIONS

1. This product may be cleaned using a mild detergent and cold water. Remove softgood from posterior stay assembly before washing. Do not use bleach or strong

- detergents. Do not wash in washing machine. machine or place in dryer.
2. Lay softgood out flat on a towel and air dry. Do not squeeze or wring out the softgood. Do not place in dryer or near heater.

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