



Instructions for Use

AIRFORM[®] Inflatable

Back Support

		4
EN	INSTRUCTIONS FOR USE	6
DE	GEBRAUCHSANWEISUNG	7
FR	NOTICE D'UTILISATION	9
ES	INSTRUCCIONES PARA EL USO	11
IT	ISTRUZIONI PER L'USO	13
DA	BRUGSANVISNING	15
SV	BRUKSANVISNING	17
FI	KÄYTTÖOHJEET	18
NL	GEBRUIKSAANWIJZING	19
PT	INSTRUÇÕES DE UTILIZAÇÃO	21
ZH	中文说明书	23

AIRFORM® INFLATABLE BACK SUPPORT

Part #	Size	Waist	U/M
209133	Small	25"-30"	1 ea
209135	Medium	30"-34"	1 ea
209137	Large	34"-38"	1 ea
209138	X Large	38"-42"	1 ea
209139	XX Large	42"-56"	1 ea

Simply add an E before the part number when ordering Hot/Cold Therapy Gel Pad.

HOT/COLD THERAPY GEL PAD

Part #	Size	U/M
GB28	Universal	1 ea

SUSPENDERS

Part #	Size	U/M
10912	Clip-On	1 ea
10913	Sewn-In	

ENGLISH

USES

Muscular low back ache, acute strains, chronic low back pain, and lumbar support

CAUTION

1. Patient should be instructed to inspect their back and this product before and after each use.
2. Do not over-inflate air bladders. No more than ten pumps should be necessary.
3. To inflate/deflate, tighten the valve until secure. Do not twist the valve too tight. Refer to sticker on knob. Release air if the brace seems too tight.

PHOTO CAPTIONS

- A. Contact closure
- B. Lumbar support
- C. Double pull contact closure strap
- D. Hot/cold therapy gel pad(optional)
- E. Air bulb

FITTING INSTRUCTIONS

1. If using hot/cold therapy, position gel pad on top of lumbar support.
2. Make sure the lumbar support is centered in the middle of the back and buttocks.
3. Fasten the contact closures to secure the support into place.
4. For additional support, pull the double pull contact closure straps from the left side and adhere in front of the abdominal area. Repeat on the right side
5. To inflate: press and release the air bulb for inflation. Repeat until the desired level of support is obtained.
6. To deflate: twist the plastic knob counterclockwise. To reinflate, twist the plastic knob clockwise until secure.

HOT/COLD THERAPY

1. For cold therapy, remove gel pad from brace and place into a freezer for approximately 20 minutes to 1 hour. To reduce swelling, apply directly to the skin for 20 minutes.
2. For hot therapy, remove gel pad from brace and place in microwave oven for 20 second intervals at low setting until hot.
Do not use oven or flame to heat the gel pad. Do not apply heated gel pad directly to skin. If too hot, allow to cool before applying.
3. Once gel pad is hot/cold, use contact closure tabs to adhere it into the brace to continue therapy. Diabetics and patients with vascular insufficiencies must always consult a physician before using hot/cold therapy.

WASHING INSTRUCTIONS

1. The brace may be cleaned using a mild detergent and cold water. Do not wash in washing machine.
2. Lay brace out flat on a towel and air dry. Do not squeeze or wring out the brace. Do not place in dryer. Do not place near heater. Do not use bleach or strong detergents.

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EN – Caution: This product has been designed and tested as a single use product, and is not intended for multiple use. In addition to increasing the risk of infection, used product components may contain small unnoticeable defects that may lead to non-optimal performance, and even mechanical failure. If any problems occur with the use of this product, immediately contact your medical professional.

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