Instructions for use

AFO Dynamic™

Life Without Limitations®

.....	3
ⓔⓃ Instructions for use.....	5
ⓓⓔ Gebrauchsanweisung.....	7
ⓕⓓ Guide de fabrication.....	9
ⓔⓢ Instrucciones para el uso.....	11
ⓔⓣ Istruzioni per l'uso.....	13
Ⓝⓞ Bruksanvisning.....	15
ⓓⓐ Brugsanvisning.....	17
ⓜⓕ Bruksanvisning.....	19
ⓕⓙ Käyttöohjeet.....	21
ⓔⓛ Οδηγίες χρήσης.....	23
Ⓝⓕ Gebruiksaanwijzing.....	25
ⓔⓣ Instruções de Utilização.....	27
ⓙⓐ 取扱説明書.....	29
ⓜⓗ 中文说明书.....	31
ⓕⓞ 사용 설명서.....	33

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



The following kits include one AFO without padding or straps

Part #	Product Name	Side	Size
O-101011	Össur AFO Dynamic	R	XS
O-101021	Össur AFO Dynamic	L	XS
O-101012	Össur AFO Dynamic	R	S
O-101022	Össur AFO Dynamic	L	S
O-101014	Össur AFO Dynamic	R	M
O-101024	Össur AFO Dynamic	L	M
O-101016	Össur AFO Dynamic	R	L
O-101026	Össur AFO Dynamic	L	L

The following kits include one AFO, padding and straps

Part #	Product Name	Side	Size
O-102011	Össur AFO Dynamic & Acc.	R	XS
O-102021	Össur AFO Dynamic & Acc.	L	XS
O-102012	Össur AFO Dynamic & Acc.	R	S
O-102022	Össur AFO Dynamic & Acc.	L	S
O-102014	Össur AFO Dynamic & Acc.	R	M
O-102024	Össur AFO Dynamic & Acc.	L	M
O-102016	Össur AFO Dynamic & Acc.	R	L
O-102026	Össur AFO Dynamic & Acc.	L	L

The following kits include only padding and straps

Part #	Product Name	Size
O-102101	Össur AFO Dynamic & Accessories	XS
O-102102	Össur AFO Dynamic & Accessories	S
O-102104	Össur AFO Dynamic & Accessories	M
O-102106	Össur AFO Dynamic & Accessories	L
O-105011	Össur AFO Dynamic & Accessories	XS
O-105012	Össur AFO Dynamic & Accessories	S
O-105014	Össur AFO Dynamic & Accessories	M
O-105016	Össur AFO Dynamic & Accessories	L

ENGLISH

PARTS LIST

- A. Proximal velcro strap
- B. Open structure tibia shell
- C. Lateral support and distal velcro strap
- D. Medial Strut
- E. Tapered footplate

INSERTION INTO SHOE

- If present, remove the insert from the shoe.
- Insert the footplate of the AFO into the shoe, if necessary grind the footplate to fit the shoeshape or adjust the shape to the contours of the foot.
- Avoid grinding excessive material from the lateral side of the footplate ensuring sufficient space between the medial strut and the medial arch of the (Drop-) foot (figure 1).
- Ensure stability of the footplate on the shoe sole throughout the entire length of the shoe.
- Make sure the footplate and (or) the insert cover the entire length of the shoe to avoid AP displacement (figure 2)
- If available, place the insert on top of the footplate.
- Have the patient don the shoe and the AFO.

FITTING THE AFO TO THE PATIENT

- Shift the AFO anterior or posterior to obtain a close fit between the distal part of the anterior shell and the lower part of the tibia (figure 3).
- If the distal section of the anterior shell is too far from the leg, it may be because the heel section of the footplate is too long. Trim the footplate at the heel end and add length to the toe part of the footplate.
- The orthosis is now ready to be fitted to the patient.

DONNING THE AFO

- The patient is in seated position with knees in 90° flexion and ankle in neutral position
- Always insert the AFO into the shoe prior to donning it to the patient.
- Align the distal section of the anterior shell to obtain close contact to the lower part of the tibia and secure this position by fastening the velcro (figure 4).
- Fixate the lower strap ensuring the position of the lower part of the shell to the tibia.
- The C(P)O is advised to determine the correct distance between the two anterior support members alongside the bony prominence of the tibia in order to ensure equal pressure distribution. Apply the strap over the 2 frontal Velcro-parts from medial to lateral (figure 5).
- Once this position is determined it doesn't require adjustment in the future. The patient can now apply the strap around the calf using one hand only (figure 6).
- Point out to the patient that doffing the orthosis, it is sufficient to relief this part of the strap, maintaining the position of the anterior support members.

CUSTOMIZING THE AFO

The following procedures are not part of the standard fitting of the orthosis and will only be carried out according to specific needs of the patient.

1. Standardized padding material is available. This can be applied according to the patient's wishes. It is advised to provide padding on the anterior shell alongside the tibia and on the medial strut covering the malleoli.
2. The use of the lower lateral support member strongly depends on the degree of voluntary control and the selection of the shoe. In case of sufficient control it is advised to remove the lateral support member according to the line shown in picture 7.
3. If the distal section of the anterior shell causes excessive pressure on the leg, the footplate should be trimmed at the toe. The lack of length at the heel is compensated by correct placement of the insert.
4. The Carbon fiber orthosis is fabricated to accommodate approximately 1,7cm (6/9 inch) heel-height. If the shoe has a significantly higher or lower heel, add or reduce the heel height by adding on material under the footplate.
5. It is recommended to use a well-fitting off-the-shelf or a custom-made insole in combination with the Dynamic AFO. This is especially the case with excessive pronators where the insole can prevent the medial malleolus and the navicular bone from getting in contact with the medial strut

of the AFO.

6. To relieve contact between the medial strut and the apex of the medial malleolus, grind the edge of the strut according to the shape shown (figure 8).
 - you may grind up to 1 cm into the strut
 - stay at least 7 cm from the footplate

GENERAL REMARKS AND ADDITIONAL ADVICE:

1. To shorten the footplate, cut or grind off the appropriate length. Be certain to buff the edges to eliminate sharp edges. To lengthen the footplate, use foot insert of choice, and fixate it to the foot plate.
2. To adjust the biomechanics of the foot, the use of the correct insert is advised.
3. **Contraindications:** Do not use the AFO in case of excessive ankle instability or in case excessive deformities to the foot where pressure relief cannot be solved by additional padding. It is contraindicated to use the AFO in case of fluctuation oedema and in case of severe spasticity.
4. The AFO cannot be reshaped by heating. Heating will cause delaminating and negatively alter the dynamics of the orthosis.
5. The right footwear is critically important for proper function of the AFO. The following criteria should be met when selecting the right footwear to use with this AFO:
 - Firm heel counter - for proper control of the rear foot
 - Firm shank - to take stress off the AFO
 - Laced or Velcro fastening - for easier donning and doffing, and to allow adjustable compressive support over the mid-foot
 - The shoes have to be wide enough to allow proper fitting of the AFO inside the shoe
 - Removable insole - to allow space for the AFO footplate with modifications
 - Approximately 6/9“ (1,7 cm) heel height - as a starting point to control the knee extension

For optimal comfort and function as well as proper customization of the ÖSSUR AFO: Dynamic, the use of an orthopaedic insole is recommended.

WARNING

With the use of this type of device, it is possible to receive various types of skin irritation and or pressure points. If you experience any unusual discomfort while using this product, consult your health care provider immediately.

This product has been designed and tested based on single patient usage and is not recommended for multiple patient use. If any problems occur with the use of this product, for example, if unusual movement or product wear is detected in a structural part of the product, immediately discontinue use of the device and contact your medical professional. Customization beyond that specified in this IFU is not recommended.

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AFO Dynamic® is covered by one or more of the following patents. Other U.S. and foreign patents pending. ® indicates trademark registration in U.S.A. and selected countries only.