ENGLISH

CONTENT: (Figure 1 & 2)

1. Valve housing
2. Rubber membrane
3. Cap
4. Cap key
5. Socket insert
6. Set screw
7. Foam cap

TOOLS REQUIRED

- 8.5mm drill bit
- M10x1.5 tap
- Hacksaw
- Caliper
- 5mm hex key
- Silicone adhesive, Össur recommends: Össur Silicone Adhesive I-400310 Philips screwdriver

VALVE INSTALLATION

Directly to the socket.

1. Fabricate the socket as usual.
2. Mark the position of the valve onto the socket.
3. Using an 8.5mm (5/16") drill bit, make a hole perpendicular to the socket wall.
4. Tap the hole using M10x1.5 tap.
5. Ensure the cap and rubber membrane have been removed prior to shortening the nozzle to prevent debris lodging in the airways.
6. Measure socket wall thickness at insertion point with a caliper, and shorten the nozzle accordingly using a hacksaw and smooth edges.

IMPORTANT: Do not shorten the nozzle beyond 4 mm (3/16") as it may affect the seal and increase the risk of the valve breaking off. It is recommended that the socket wall thickness at valve insertion point is to be equal or more than 4mm (3/16") in thickness. (Figure 3)

Using the socket insert. (Figure 2)

1. The socket insert is used as a fixation device for the valve, and should always be used when fabricating a flexible inner socket.
2. Prepare the plaster cast as usual.
3. Fix the socket insert to the cast using the set screw. Ensure it sits flush to the cast.
4. Plaster may have to be removed or added.
5. Apply the foam cap to the set screw.
6. Fabricate the socket as usual.
7. Grind material away over the socket insert until foam cap is fully exposed.
8. Unscrew the set screw to remove socket from cast.
9. Ensure the cap and rubber membrane have been removed prior to shortening the nozzle to prevent debris lodging in the airways.

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- 8,5-mm-Bohrerspitze
- M10x1,5 Gewindeschneider
- Metallsäge

10. Measure socket wall and socket insert thickness at insertion point with a caliper, and shorten the nozzle accordingly using a hacksaw and smooth edges.
- IMPORTANT:** Do not shorten the nozzle beyond 4 mm (3/16") as it may affect the seal and increase the risk of the valve breaking off. (Figure 3)

ASSEMBLY.

1. Apply silicone adhesive to the threads and housing, where it comes in contact with the socket.
- NOTE:** It is not recommended to use epoxy or similar as it can compromise the seal.
2. Install valve and hand tighten with a 5mm hex key.
3. After installation distribute silicone adhesive around the valve / socket contact point to ensure an air tight seal.
4. Ensure the valve parts and airways are free from debris before final assembly.
5. Replace rubber membrane to cap. Ensure the concave surface is facing the cap cavity. (Figure 4)
- NOTE:** When reassembling ensure the membrane is ON the cap and NOT in the housing.
6. Screw the cap on using the cap key.

DISASSEMBLY AND CLEANING

1. Unscrew the cap using the cap key.
2. Remove the rubber membrane carefully.
3. Clean all parts with water and mild detergent.
4. Hold down the release button while cleaning airways.

REASSEMBLY

1. Replace rubber membrane to cap. Ensure the concave surface is facing the cap cavity. (Figure 4)
2. Screw the cap on using the cap key.
- NOTE:** When reassembling ensure the membrane is ON the cap and NOT in the housing.

PRODUCT WARRANTY

ÖSSUR offers 6 months limited warranty on manufacturing faults.

A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

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EN – Caution: Ossur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Ossur products are used with other recommended Ossur components. If an usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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©Copyright Össur 2017
Össur Deutschland GmbH