Instructions for use

UltraWrap

Life Without Limitations®

.....	3
ⓔⓃ Instructions for use.....	5
ⓓⓔ Gebrauchsanweisung.....	6
ⓕⓗ Guide de fabrication.....	7
ⓔⓈ Instrucciones para el uso.....	8
ⓔⓣ Istruzioni per l'uso.....	9
Ⓝⓞ Bruksanvisning.....	10
ⓓⓐ Brugsanvisning.....	11
Ⓢⓕ Bruksanvisning.....	12
ⓕⓙ Käyttöohjeet.....	13
Ⓝⓕ Gebruiksaanwijzing.....	14
ⓓⓣ Instruções de Utilização.....	15

EN - Caution: This product has been designed and tested based on single patient usage and is not recommended for multiple patient use. If any problems occur with the use of this product, immediately contact your medical professional.

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




ENGLISH

FITTING INSTRUCTIONS

1. Place opened garment behind back of knee.
Note: garment should be stretched lightly before first application.
2. Wrap the bottom portion of the garment around the calf and secure loosely.
3. Wrap the top portion of the garment around the thigh and secure loosely.
4. Position sleeve so the opening is centered over the patella. Adjust each of the 4 contouring pull points until the garment fits snugly and comfortably.
5. Fasten the lower strap snugly across the calf.
6. Fasten the upper strap snugly across the thigh.

CHANGING EXTENSION STOPS

7. Remove hinges from garment. Using a short Phillips screwdriver, remove the screw from the center of the hinge plate .
8. Locate extension stop positioned in the hinge between the upper and lower arms. With arms slightly flexed, push the extension stop from the back of the hinge through to the front, causing the stop to pop out. Feed the hole end of the desired stop through the hinge from the front to the back, with the curved side facing down and the flat edge facing up.
9. Make sure that the sharp extrusions on the flat end of the stop are pointed downward. The flat end of the stop should be visible from the front. Align the hole in the center of the hinge plate with the hole in the extension stop. Replace and tighten the screw through the hinge plate and extension stop.

WASHING INSTRUCTIONS

- Remove hinges and close hook/loop fasteners before washing.
- Use warm, soapy water, rinse thoroughly and line dry only.
- Do not machine wash or dry.
- For best performance results, rinse after every use.

WARNING: this product is designed for warmth and compression only and will not prevent injury. If you experience skin irritation, pain, or discomfort while using this product, remove immediately and consult your health care provider

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