Instructions for use

Rebound™ Ankle Brace

Life Without Limitations®

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



READ THE FOLLOWING INSTRUCTIONS CAREFULLY BEFORE USING THE DEVICE

Caution: This product has been designed and tested based on single patient usage and is not recommended for multiple patient use.

- Ensure patient performs circulation checks. If a loss of circulation is felt, or if the brace feels too tight, adjust strap tension to reduce compression to a comfortable level. If discomfort continues, patient should contact their medical professional immediately.
- If patient experiences discomfort, for example pressure points, adjust strap tension. If discomfort remains, patient should contact their medical professional for adjustments.

INDICATIONS AND USES

Ankle sprain grade I, II, and III, chronic ankle instability, high ankle sprain (syndesmosis), transition out of a walker and/or cast after fracture or surgery, posterior tibialis dysfunction, and prophylactic use.

CONTRAINDICATIONS

Ankle fracture, excessive swelling, immobilization of dorsiflexion and/or plantarflexion

APPLICATION

1. Insert product into shoe. Unfasten and open dual closure straps (A) (Fig. 1)
2. Insert foot into brace and shoe. Tie shoe laces.
3. Cross, wrap, and fasten dual closure straps (A) (Fig. 2)
4. Cross the stability strap (C) across the top of the foot (Fig. 4)
5. Snugly pull the stability strap (C) around and behind the ankle (Fig. 5)
6. Fasten the rubber strap end (D) down on the stability strap (C) (Fig. 6)

Note:

- This product is intended to be worn with an athletic sock.
- For quick re-application, leave in shoe.
- The stability strap (C) is able to be trimmed with scissors but may need to be a different length if used in a different stability configuration. After trimming, the rubber strap (D) must be re-attached.

ADJUSTMENT FOR WIDTH

1. Un-snap and remove the width adjustment snaps (B)
2. Place into slots “1” for small calves, “2” for medium calves (standard), “3” for large calves (Fig. 3)
3. To reduce bulk, trim the liner at the appropriate indentation under the dual closure straps (A).

Note: If you have purchased the width extender accessory kit, place into slots “1” for X-Large calves, “2” for XX-Large calves, “3” for XXX- Large calves. (Extender for XL-XXXL calves available for additional charge)

INSTRUCTIONS FOR HEAT MOLDABLE FOOTPLATE

1. Remove the soft-good.
2. Using a heat gun, apply approximately 350° F (175° C) of heat evenly to the outside surface (surface not in contact with the sock) of the desired area of the footplate. Keep the heat gun 6” (15 cm) from the surface. Apply heat for approximately 20-30 seconds until the plastic is malleable but not melted. Avoid heating rivets – may burn skin.

- Using gloves, manually form the footplate into the desired form. Allow sufficient cooling time prior to applying device. Do not attempt to form components against bare skin.

CLEANING INSTRUCTIONS

Hand wash soft-good in warm water with mild soap. Air-dry only. Wipe plastic clean with a damp towel and air dry. Do not machine wash/dry.

OPTIONAL ACCESSORIES

To order width extender or stability strap kit, contact your healthcare provider or Össur customer service

STABILITY STRAPPING CONFIGURATIONS

Note: Square hook sticker must be applied to the bottom of the footplate for all configurations.

- EZ Stabilizer:** (default configuration): The strap initiates from bottom of the footplate (Fig. S1), wraps up over the inside of the foot, and then around the ankle (Fig. 6).
- EZ Stabilizer Light:** The strap initiates from bottom of the footplate (Fig. S2), wraps up over the outside of the foot, and then around the ankle (Fig. S3).
- EZ Stabilizer X:** An additional hook sticker must be applied to the inner strut (Fig. S5). The middle of the strap initiates on the bottom of the footplate (Fig. S4), both ends cross over the top of the foot (Fig. S5), and wrap around the ankle (Fig. S6).
- Stabilizer:** An additional hook sticker must be applied to the inner strut (Fig. S7). The strap initiates from inside the inner strut (Fig. S7), wraps around the outside of the foot, under the footplate (Fig. S8), up across the top of the foot, and then around the ankle (Fig. S9). For this configuration, the brace cannot remain in the shoe for quick reapplication.
- Stabilizer Plus:** An additional hook sticker must be applied to the outer strut (Fig. S10). The strap initiates from inside the outer strut (Fig. S10), wraps around the outside of the foot, under the footplate (Fig. S11), up across the top of the foot, and then around the ankle (Fig. S12). For this configuration, the brace cannot remain in the shoe for quick reapplication.

	EZ Stabilizer Light	EZ Stabilizer	EZ Stabilizer X	Stabilizer	Stabilizer Plus
Ease of Application	+++	+++	++	+	+
Foot Stability	+	++	+++	+++	++++
Ankle ML Stability	++	++	++	++	+++

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