Instructions for use

Techform[®] Premium Casting Tape

Life Without Limitations[®]

EN - Caution: This product has been designed and tested based on single patient usage and is not recommended for multiple patient use. If any problems occur with the use of this product, immediately contact your medical professional.

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® indicates trademark registration in U.S.A. and selected countries only.

ENGLISH

Please read these instructions before use. The healthcare professional applying this casting material is responsible for the proper care of the patient. Only responsible, qualified, and experienced healthcare professionals should apply casting tape to patients.

INDICATIONS FOR USE:

TECHFORM® casting materials are to be used for casting, splinting, and external immobilization.

MATERIAL DESCRIPTION:

TECHFORM® casting materials consist of a knitted synthetic substrate impregnated with a water activated polyurethane resin.

PRODUCT STORAGE:

This product should be stored at room temperature. Storage at high temperatures will result in a reduced shelf life and tacky handling characteristics. In addition, care should be taken to avoid puncturing or tearing the pouch in storage as this will cause premature hardening in the package.

FLAMMABILITY:

Uncured polyurethane tape will burn over an open flame. Finished casts using TECHFORM® cast tape will not support combustion under normal conditions as defined by the NFPA Standard #702-1975 and the Flammable Fabrics Act (CS-191-53).

PRE-APPLICATION PRECAUTIONS:

- The healthcare professional applying this casting material is responsible for the proper care of the patient. Only responsible, qualified, and experienced healthcare practitioners should apply casting material to patients.
- The polyurethane resin will adhere to unprotected skin or clothing. Apply gloves before handling the casting material.
- Exposure of the resin to ambient air or water will begin the activation process and the tape will become rigid. Therefore, it is important that the cast tape pouch remain free from punctures until it is time to submerge the cast tape in water and apply the cast.
- As the material begins to set, it releases heat (exotherm).

The following are conditions that can increase heat to uncomfortable/ injurious levels:

1. The more layers applied to the cast, the higher the exotherm, therefore, care should be taken when applying multiple layers of cast tape.
 2. Activating the cast tape by submerging the roll in water over 80° Fahrenheit (27° Celsius) may result in a higher exotherm temperature.
 3. Applying any coverings such as wraps, dressings, and/or padding over a setting cast, or splint, or lying the cast on any surface such as a pillow during the setting process should be avoided. In addition, proper care should be exercised with children and patients with sensitive skin.
- Do not insert any materials in between the cast and the patient or padding.
 - Inform the patient that they should not get the cast wet. Skin maceration and irritation may result from prolonged and/or frequent wetting of the cast.

CAST APPLICATION:

- Immerse the cast tape roll in room temperature water. (68° F (20° C) to 75° F (24° C)) and firmly squeeze the roll three times while it is submerged in the water. Only open one pouch at a time.
- Traditional casting techniques are to be exercised when applying the casting material. 3-6 layers of cast tape are needed for proper immobilization. Weight bearing casts should consist of a minimum of four layers. (See heat precautions)
- The product should be rubbed well during application for ideal lamination of the cast tape layers. Set times vary between 3 to 5 minutes when the material is activated with room temperature water.
- Allow 30 minutes before allowing the patient to bear weight or loads on the cast.

Removal, alteration, uni/bivalving, and repairing should be instituted with care. The applicator is responsible for any alteration and/or removal of the cast or splint.

CAST REMOVAL:

Use a standard cast saw for cast removal. Care should be taken to avoid bony prominences, as well as, excessive heat build up from the cast saw blades. Avoid exposed edges of the cast to prevent scratches.

PATIENT INSTRUCTIONS:

Instruct the patient of the following information to help insure patient comfort and safety, as well as, care for the cast integrity.

- Patients should be notified by their healthcare professional as to proper care for their injury, cast instructions, mobility instructions, and how to notify their healthcare professional should problems with the cast occur.
- Patients should be encouraged to not get the cast wet. The cast may be cleaned with a damp cloth if necessary. Should the cast become wet, use a towel to remove excess moisture and then use a blower dryer (such as a hair dryer on a cool setting) to sufficiently dry the undercast materials.
- Skin maceration and irritation may result from frequent wetting of the cast without proper drying.
- Do not insert any materials into the cast. Avoid activities that may result in dirt, sand, or other foreign articles getting under the cast.
- If necessary, rough edges of the cast may be reduced with a nail file.
- The patient should immediately contact their healthcare professional if they experience; painful pressure or rubbing from the cast, their toes or fingers become numb or bluish in color, or if the cast feels too tight. In addition, the patient should notify their healthcare professional if the cast softens or breaks at any time.

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