

Quick Reference Card

CHARGING AND POWER

1. Charge the RHEO KNEE® 3 for 3 hours to ensure full charge
2. RHEO KNEE 3 can be charged ON or OFF

Important: During assembly maintain minimum distance of 3 mm (1/8") between the RHEO KNEE 3 and the socket at maximum flexion.

ALIGNMENT

1. The alignment reference line should pass through a sagittal plane bisection of the socket, the knee center and the posterior 1/3 of the Flex-foot.
2. **Note:** If stability modifications are desired after patient assessment; the reference line can be shifted up to 5 mm anterior or 5 mm posterior to the knee center.



CONNECTING

1. Ensure the RHEO KNEE 3 is switched on
2. Connect to the RHEO KNEE 3 via Bluetooth and note the COM port number
3. Open RHEOLOGIC and connect to the RHEO KNEE 3 using the corresponding COM port number

AUTO ADJUSTMENT

1. Select Auto Adjustment Tab
2. Press "ON"
3. Have the user walk 30 steps
4. Press "OFF"



EXTENSION LOCK

