ENGLISH

SETTING UP THE PCL SUPPORT SYSTEM

The PCL kit is designed to counteract against a posterior subluxation of the tibia. The thigh plate and strap apply a posterior force against the anterior thigh, while the lower straps and calf plate apply an anterior force to the tibia. To place the PCL system in the correct position on the leg:

1. Remove the upper cruciate d-rings.
2. Replace the female insert, original d-ring, and washer on the outer surface of the frame.
3. Insert the male insert (double dots on the head of the male insert) into the circular recess on the PCL d-ring.
4. Fasten down the assembly with a screw.
5. If the double dot male insert with the 2 dots is too long, then replace with the single dot male insert.
6. Repeat for the remaining upper cruciate d-ring.
7. Apply the brace.
8. Loosely apply the thigh strap with the anterior thigh plate into the PCL d-rings (straps may need to be trimmed to size).
9. Attach PCL cuff to the lower straps and center over the calf muscle. Tighten straps so that the cuff fits comfortably but snugly against the back of the leg.
10. Tighten anterior thigh strap with the leg in full extension so that the plate and pad are comfortably snug.

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©Copyright Össur 2014