
Instructions for Use

---



FLEX-RUN<sup>™</sup> Junior

with Nike Sole





Product number: FSX007x

		3
EN	INSTRUCTIONS FOR USE	4
DE	GEBRAUCHSANWEISUNG	6
FR	NOTICE D'UTILISATION	8
ES	INSTRUCCIONES PARA EL USO	10
IT	ISTRUZIONI PER L'USO	12
DA	BRUGSANVISNING	14
SV	BRUKSANVISNING	16
NL	GEBRUIKSAANWIJZING	18
PT	INSTRUÇÕES DE UTILIZAÇÃO	20
PL	INSTRUKCJA UŻYTKOWANIA	22
JA	取扱説明書	24

**BENCH ALIGNMENT**

# ENGLISH

---

## FLEX-RUN™ SYSTEM (Figure 1):

- Custom foot for recreational sports activities
- Preferred long distance running foot
- Build Height: 261mm (10 1/4"), clearance: 208mm (8 3/16")
- Maximum weight 130 kg (285 lbs)
- Endoskeletal attachment options

## CONNECTOR OPTIONS:

- Male Pyramid Foot Adapter Set CM131138
- Pyramid Adapter FWX00030 (1–6)
- Female Pyramid CM 150161 (1–6)
- Attachment plate FWX40001 must be added when using male or female adapter.
- Roughen top surface of Flex-Run (Figure 2) and attachment plate FWX40001
- Mix Siegelharz and apply to both surfaces
- Bond both parts together and clamp securely using a greased center bolt and the connector part
- Disassemble after resin has set and trim off excess
- Clean parts, assemble and torque using Loctite 410 to following specifications:

<b>Center bolt</b>	M10	41 Nm (30 lbs)
	M12	43 Nm (32 lbs)
<b>Pyramid angulation</b>	M 8	15 Nm (9 lbs)

## BENCH ALIGNMENT (Figure 3 - 4)

### SAGITTAL PLANE:

- Place the ground contact point of the foot module 25–40 mm (1–1½") in line with load line.
- Determine the appropriate socket flexion according to amputee requirements.

### CORONAL PLANE:

- Introduce appropriate socket abduction/adduction.
- Allow for up to 5° of external rotation.
- Achieve proper height approximately 25–50 mm (1–2") taller than the sound side with shoe to compensate for vertical compression.  
Note: The recommendations for bench alignment are used as a starting point. Definitive alignment will be achieved during dynamic assessment.

## DYNAMIC ALIGNMENT

Physical ability and length of residual limb will affect alignment. Best performance will be achieved with best achievable symmetry of gait. Each athlete will have her/his own optimal alignment, which depends on the individual's running style. Optimal alignment is critical in order to reduce energy expenditure and eliminate excessive muscle and ligament strain while running:

- Consult with your Össur clinician before using the running foot for the first time.
- Check for best symmetry of gait.

- Stride length can be affected by A-P alignment and vertical deflection socket to foot.
- Vertical deflection is between 25 mm (1") and 50 mm (2") on average.

### **SOLE ATTACHMENT (Figure 5 - 6)**

- Slide the Nike sole on to the foot and make sure all the clips go over the foot.
- Pull the rubber strap of the sole over the medallion and slide it under. When detaching the sole, start by pulling the strap away from the medallion, then slide the sole off the foot.
- For tighter fit of the sole or to reduce noise bond insert inside sole using Loctite, keep under pressure for 10 minutes cure

### **TROUBLESHOOTING**

Make sure socket fit and suspension are ideal. Symmetry in running should be assessed in the optimal individual speed phase.

***If the duration of loading the prosthesis is longer than the contralateral side:***

- Foot module is too soft
- Prosthesis is too short
- Foot is too far anteriorly aligned

***Duration of loading the prosthesis is shorter:***

- Foot module is too stiff
- Prosthesis is too long
- Foot is too far posterior aligned

### **TRANSFEMORAL ATHLETE CONSIDERATIONS**

Moving the foot in A-P position will affect hip flexion, stride length and trunk motion.

### **WATERPROOF**

A waterproof device allows use in a wet and/or humid environment and temporary (<30 mins.) submersion of up to 1m in fresh water. Exposure to fresh water under these circumstances shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

The components used with the foot need to be waterproof as well.

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**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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