








Instructions for Use

LP VARI-FLEX
LAMINATION ADAPTER

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	10
ES	Instrucciones para el uso	12
IT	Istruzioni per l'uso	15
DA	Brugsanvisning	18
SV	Bruksanvisning	20
NL	Gebruiksaanwijzing	22
PT	Instruções de Utilização	25
JA	取扱説明書	28

ENGLISH

LP Vari-Flex Lamination Adapter Kit is exclusively to be used with the LP Vari-Flex foot module for lower limb prosthetics. The VLX00001-3 LP Vari-Flex Lamination Adapter Kit is available for users using LP Vari-Flex in indicated sizes and categories:

1. size 22–24, category 1–5
2. size 25–30, category 1–9

The VLX00001-3 LP Vari-Flex Lamination Adapter Kit is rated up to 166 Kg (365 lbs) – category 9.

EACH KIT INCLUDES: (Figure 1)

- 1× lamination adapter,
- 4× attachment bolts,
- 1× strip of sealant tape,
- 1× tube of Loctite 410

BENCH ALIGNMENT (Figure 2+3)

- Fit foot with selected cover.
- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/extension and abduction/adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram (the mark on the brim of the foot cover matches the one third portion).

FABRICATION PROCEDURES

- Fill the space between the socket and the foot with a material of choice such as potting epoxy, PU foam or wood, depending on weight and impact level of user.
- Bond the VLX30201-3 lamination adapter to the filling material.
- Wrap synthetic casting tape or similar tape and proceed with dynamic alignment.
- When forming the filling material, avoid sharp edges and horizontal sections as that reduces the strength of the socket (**Figure 8+9**).

DYNAMIC ALIGNMENT

The heel stores energy after initial contact, slowly releasing it at midstance. The forward momentum generated by this action results in the toe being loaded for optimum energy release at terminal stance. The heel to toe action can be influenced by

- A–P positioning of foot
- Dorsi-Plantarflexion
- Heel Stiffness

LAMINATION TECHNIQUE

- Invert socket and secure in place.
- Clean the groove of the aluminum adapter with Acetone or similar to ensure maximum bonding of the fibers to the aluminum.
- If foam or another porous material has been used as the filling material it is
- recommended to coat it with resin to ensure that resin will not be drawn into the foam when laminating the socket.

- Position the appropriate amount of carbon braids or chosen reinforcement. For better attachment it is recommended that the braids are layered over the aluminum piece.
- Tie off the overhanging braids around the base groove using a heavy thread. **ENSURE THAT THE TIE STAYS TIGHT (Figure 5).**
- Fold the overhanging braids over the socket.
- Finish with a stockinette for cosmetic finish, if desired.
- Seal holes thoroughly using the sealant tape provided and apply grease on the bottom attachment surface of the aluminum adapter (alternatively seal the surface with PVC tape). **(Figure 4)**
- Pull a PVA bag and pour in resin.
- Tie off excess resin before it reaches the gel state. Allow resin to harden.
- Grind off the excess cured resin around the edge of the aluminum adapter.
- Expose the attachment surface. Care must be taken to not grind the aluminum piece itself!
- Assemble and torque attachment bolts according to specifications using Loctite 410. To avoid stripping the hex-slot of the bolt, make sure to use sharp hexkeys and ensure full engagement of key to hex slot of the bolts. See also instructions for LP Vari-Flex **(Figure 6+7).**

IMPORTANT

It is important that the torque specifications are met.

NOTE

Do not reuse adapter attachment bolts, which have been torqued according to specifications or if Loctite has been applied to them. New bolts must be used.

Table 1: Suggested reinforcement

Size and category	Layers of carbon fiber
Size 22–24, cat 1–5	Minimum 4*
Size 25–30, cat 1–6	Minimum 6*
Size 25–30, cat 7–9	Minimum 6*

* The appropriate number of layers depends on user activity and weight and must be evaluated by the CPO. It is recommended to use heavy carbon fiber braids.

Table 2: Torque specifications (Figure 6+7)

Size and category	Adapter attachment bolts (Nm)	Heel plate attachment bolts (Nm)
Size 22–24, cat 1–5	22	14
Size 25–30, cat 1–6	34	14
Size 25–30, cat 7–9	34	28

A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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