
Instruction for Use

BALANCE KNEE OM8

		5
EN	Instructions for Use	7
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	16
IT	Istruzioni per l'uso	19
DA	Brugsanvisning	22
NL	Gebruiksaanwijzing	25
PT	Instruções de Utilização	28
RU	Инструкция по использованию	31

ALIGNMENT

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


Loodlijn

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BALANCE KNEE OM8

The Balance Knee OM8 is a polycentric knee joint with a stable aluminium frame construction. The joint is suitable for users with low functional levels, with a high to very high need for safety, and a maximum body weight of 136 kg.

INDICATIONS FOR USE

The Balance Knee OM8 is designed for patients with transfemoral and knee disarticulation amputation.

CONTRAINDICATIONS FOR USE

N/A.

TECHNICAL DATA

Mobility class	1, 2
Max. body weight	136 kg
Material	AL
Amputation level	Above-knee
Overall fitted height	approx. 128 mm
Effective fitted height	approx. 76 mm
Fitted height	approx. 12 mm
Weight	approx. 445 g
Prox. connection	Pyramid
Ground clearance	approx. 1.6 cm
Dist. connection	Integrated 30 mm tube receiver
Swing phase	Integrated extension spring, friction adjustment
Stance phase	Axle configuration, fulcrum configuration
Axes	Polycentric
Knee flexion without socket	approx. 160°

ALIGNMENT

Stance phase stability

Axial alignment can be changed using the 5 mm Allen screw (25). The joint is locked when delivered (**fig. 1**).

- Turn the 5 mm Allen screw (25) clockwise to make the alignment more dynamic.
- Turn the 5 mm Allen screw (25) anticlockwise to make the alignment safer.

CAUTION: make small adjustments and check the result! (fig. 2)!

Swing phase control

Adjusting the extension spring (**fig. 3**):

- Turn the adjusting screw (1) clockwise with the help of a screwdriver: the extension function increases.
- Turn the adjusting screw (1) anticlockwise with the help of a screwdriver: the extension function decreases.

Adjusting the axle friction (**fig. 4**):

- Turn the adjusting screw (15) with the help of a 5 mm T spanner clockwise, axle friction increases.
- Turn the adjusting screw (15) with the help of a 5 mm T spanner anticlockwise, axle friction decreases

SAFETY WARNINGS

- Adapters used on the proximal connection must always be tightened with the set torque of the relevant fit component. Use a torque wrench.
- Adapters (e.g. tube adapters) used on the distal connection must be inserted down to the end stop in the tube receiver of the knee joint.
- Tighten the screw (6) on the tube receiver (5) to a torque of 16 Nm and secure with Loctite. Use a torque wrench (5 mm Allen).
- Observe the alignment guidelines and adjustment instructions.
- Do not use any screws apart from those described above.
- Please do not use any talcum since this could block the joints, above all in connection with lubricants.
- Do not use compressed-air for cleaning. We recommend the use of a brush.

Observe all the above points, otherwise the guarantee becomes null and void!

MAINTENANCE

The joint must be inspected and repaired if necessary every 6 months!

Inspect for:

- the alignment
- the screw connections
- the suitability of the patient (e.g. max. body weight, mobility)
- max. weight and mobility
- loss of lubricant
- damage to the joint and accessories
- soiling of the bushes
- joint play A/P and/or M/L
- contamination with talc
- the extension effect of the extension spring (2)
- the extension stop (23)

CARE

- Clean the joint with a soft cloth moistened with a little mild benzene. Do not use any more aggressive cleaning agents because these can damage the seals and bushes.
- Do not use compressed air for cleaning! Compressed air can force dirt into the seals and bushes. This can lead to premature damage and wear and tear.

LIABILITY

In the event of a claim: Össur will only consider complaints submitted with a copy of the delivery note or invoice from Össur and exact details of the reason for returning the joint. A manufacturer can only be held liable for the failure of his own fittings. Any liability of the manufacturer in excess of this is only possible if his fittings are demonstrably the cause of the damage or loss of function of fittings from other manufacturers.

CE MARK, COMPATIBILITY

All modular parts have been tested to EC Directive 93/42/EEC and bear the CE mark.

General Conditions:

This product is recommended for use by one patient only. Product liability pursuant to the Medical Devices Act becomes null and void if it is used for treating more than one patient.

The compatibility of Össur products with fittings from other manufacturers is guaranteed as long as the following points are satisfied:

- The fittings are only used in compliance with their intended purpose (see Instructions for Use).
- If fittings with different max. weight limits are used, the maximum loading corresponds to that of the weakest fitting.

- The use of tested individual components does not relieve the prosthetist from his duty to check by all possible means the combination of fittings for their suitability, correct assembly and safety.
- If there are any indications that a certain combination of fittings does not guarantee the required level of safety, then the fittings must not be combined.
- The prostheses must be built in accordance with the generally recognised professional rules of the orthopaedic technical trade.
- If the patient's situation changes (e.g. body weight, degree of activity, etc.) the prosthesis must be checked immediately.
- If the prosthesis has been subjected to disproportionately high stresses, e.g. during a fall, it must be inspected for any possible damage immediately by a specialist workshop.
- Safety-relevant regulations for individual fittings (e.g. maintenance intervals, see Instructions for Use) must be observed.

Combinations of different fittings:

The above general conditions apply, in addition:

- Use only fittings that comply with the requirements of ISO 10328.
- Use only fittings that comply with the requirements of the Medical Devices Act.
- The fittings must be designed for use with standard pyramid or 30 mm systems, and must be compatible with each other.

USEFUL LIFE

The applicable provisions of the latest Össur Product Catalogue apply. Maximum useful life of the Balance Knee OM8 is 36 months.

REPLACEMENT PARTS/ACCESSORIES

Only the listed replacement parts and accessories are available.

No.	Description	Article no.
1	Adjusting screw	1.980.100
2	Spring	1.980.101
6	Clamp screw	1.980.104
15	Flexion screw	1.980.106
23	Extension stop	1.980.107
24	Washer	1.980.108
25	Thread Screw	1.980.109

To have a joint inspected please send the joint for repair to Össur. We will provide you with a service joint for the duration of the repair.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Össur products are used with other recommended Össur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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