
Instruction for Use

OP2 KNEE

		5
EN	Instructions for Use	7
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	22
DA	Brugsanvisning	26
NL	Gebruiksaanwijzing	29
PT	Instruções de Utilização	32

Aufbau

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







OP2 KNEE

Art.Nr. 1.724.120

The OP2 Knee is a polycentric knee joint with a sturdy aluminium frame construction. It is suitable for users in Mobility classes 2 and 3 and is approved up to a body weight of 125 kg. The OP2 Knee is for users with normal demands on walking (steady, low walking speed) and an increased safety requirement.

INDICATIONS FOR USE

The OP2 Knee is designed for patients with transfemoral amputation.

CONTRAINDICATIONS FOR USE

N/A.

TECHNICAL DATA

Material	aluminium
Amputation height	thigh
Weight limit	125 kg
Tot. installation height (1)	approx. 190 mm
Effective overall height (2)	approx. 137 mm
Overall height (3)	approx. 12 mm
Weight	approx. 650 g
Ground Clearance	approx. 16 mm
Proximal connection	pyramid
Distal connection	integrated 30mm pylon
Swing phase	pneumatic
Stand phase	adjustable axis geometry
Axes	polycentric
Knee bend angle without shaft	approx. 160°

ADJUSTMENT

The OP2 Knee is delivered ready for use in the basic position ex works.

Start testing any adjustment of the joint slowly with the user during the dynamic walking test and check whether the basic position suits the needs of your user. If it does not proceed in the manner described below.

Safeguarding the stand phase

The axial arrangement can be altered with the 5 mm Allen screw (29) (**Fig. 1**). The joint has a safe setting on delivery.

- Turn the 5 mm Allen screw (29) clockwise for a more dynamic setting.
- Turn the 5 mm Allen screw (29) anti-clockwise for a safer setting.

Caution: Carry out all adjustments in small steps and check the results (Fig. 2-4)!

Controlling the swing phase

The pneumatic system is delivered in its basic position. In the basic position each of the two valves is open by two revolutions. You can carry out adjustments to this basic position to match the individual gait pattern. First adjust the flexion to achieve an acceptable walking pattern. You can then harmonise the gait pattern by means of the extension.

Flexion too much

Turn the flexion valve clockwise with a 3 mm Allen key. (Fig. 7)

Flexion not enough

Turn the flexion valve anti-clockwise with a 3 mm Allen key. (Fig. 8)

Bumper too hard

Turn the extension valve clockwise with a 3 mm Allen key. (Fig. 9)

Extension too slow

Turn the extension valve anti-clockwise with a 3 mm Allen key. (Fig. 10)

Caution: The extension may only be dampened enough so that a full extension is always possible!

SAFETY INSTRUCTIONS

- The adapter used (e.g. tube adapter) on the lower connection must be slid up to the stop in the knee joint pylon clamp.
- The pylon clamp (8) screw (5) must be tightened with a torque of 16 Nm and secured with Loctite. Use a torque wrench (5 mm hexagon socket).
- The assembly and adjustment instructions must be observed.
- No screws other than those described here may be opened.
- The use of talcum powder may block the joint, especially in connection with lubricants. We thus ask that you refrain from using talcum powder!
- Do not use compressed air for cleaning. We recommend the use of a brush.

All of the points mentioned above must be observed since the warranty will otherwise be forfeited.

SERVICING

The joint should be serviced at least every 6 months.

The following should be checked during the service:

- the assembly
- the screwed connections
- the patient care requirements (e.g. weight limit, mobility)
- lubricant leaks
- damage to the joint and accessories
- dirt on bushings
- AP and/or ML play in the joint
- soiling due to talcum powder
- the extension effect of the extension assist (4)
- the bumper (31)

CARE

- Clean the joint with a soft cloth steeped in a little mild petroleum ether. Do not use more aggressive cleaner since this could damage the seals and bushings.
- Do not use compressed air for cleaning! Compressed air can bring dirt into seals and bushings which may then lead to damage and premature wear.

LIABILITY

The following applies in the event of damage or loss: Össur will only process the complaint with a copy of the delivery note or invoice from Össur with detailed reasons for why the joint is being returned.

Each manufacturer can only be held liable for the failure of their own fit components. Any further liability on the part of the manufacturer is only possible if it can be proven that his fit components were the cause of the damage or malfunction of fit components from other manufacturers.

CE SYMBOL/COMPATIBILITY

All modular parts have been tested in accordance with EC Directive 93/42/EEC and bear the CE symbol.

General terms and conditions

This product is only intended for use with one patient. If it is used to treat more than one patient the product liability in the intendment of the German medical products act will be forfeited.

Össur products are guaranteed to be compatible with other fit components if the following points are met:

- The fit components are used only for their intended purpose (pay attention to the relevant instructions for use).
- If fit components with different max. weight limits are used, the max. load of the weakest fit component applies.
- The use of tested individual components with a CE symbol does not release the technician from his obligation to check the combination of fit components within the scope of his possibilities with respect to their suitability, correct assembly and safety.
- If there are any indications that a combination of fit components does not offer the necessary safety then these fit components may not be combined.
- The prosthesis must be constructed in accordance with the generally accepted rules of the orthopaedic technology trade.
- If the patient prerequisites change (e.g. body weight, mobility, etc.) the appliance must be checked immediately.
- If the appliance has been exposed to disproportionately high stress (e.g. fall) it has to be examined immediately by a specialist workshop for any possible damage.
- Any regulations relating to individual fit components that are relevant for safety (see the relevant instructions for use) must be observed.

Combination of different fit components

- In addition to the general terms and conditions listed above, the following also applies:
- only use fit components that satisfy the requirements of ISO 10328
- only use fit components that satisfy the requirements the MPGs
- the fit components must be designed for the common pyramid or 30 mm system and be matching.

USEFUL LIFE

The maximum useful life of the OP2 Knee is 36 months.

SPARE PARTS, ACCESSORIES

The following spare parts are available:

No.	Material	Article number
4	Extension assist spring	1.970.252
5	Pylon clamp screw	1.970.866
6	Extension adjustment valve needle	1.980.205
7	Flexion adjustment valve needle	1.980.206
29	Bumper adjustment screw	1.980.109
30	Metal washer	1.980.108
31	Bumper	1.980.107

Only the listed articles are available as spare parts and accessories.

Please send the joint to be repaired to Össur for examination. We will provide a service joint for the duration of the repair.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Össur products are used with other recommended Össur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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