

Instructions for Use

DIRECT SOCKET TOOL KIT


		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	12
IT	Istruzioni per l'uso	14
NO	Bruksanvisning	17
DA	Brugsanvisning	19
SV	Bruksanvisning	21
EL	Οδηγίες Χρήσης	23
FI	Käyttöohjeet	26
NL	Gebruiksaanwijzing	28
PT	Instruções de Utilização	30
PL	Instrukcja użytkowania	33
CS	Návod k použití	36
TR	Kullanım Talimatları	38
RU	Инструкция по использованию	40
JA	取扱説明書	43
ZH	中文说明书	45
KO	사용 설명서	47

3

Direct Socket pressure casting	New amputee	Fleshy	Muscular	Bony
mmHg	40-50	60	70	80

ENGLISH

PRODUCT DESCRIPTION

Direct Socket Tool Kit is a single chamber pressure casting system providing pressure to shape soft tissue. The single chamber pressure system is designed to optimally distribute pressure. The chamber is reinforced with a matrix for durability, and to control elongation.

Direct Socket Tool Kit provides an effective and clinically validated option for socket provision. It provides the prosthetist with tools to repeatedly produce a reliable and well-fitting socket using the Direct Socket process. Optimal axial and transverse stabilization can be obtained by using the Iceross Casting Liner (I-09DCXX) in combination with the Direct Socket Tool Kit.

Please refer to the Instructions for Use of the Direct Socket System for handling of Resin Injection Tool, Distal Attachment, Taping Ring and Molding Tube.

Direct Socket Tool Kit contents:

- Icecast Pump (D-040510)
- Icecast Bladder (D-040520)
- Resin Injection Tool (M-400201(200 ml) & 1950504 (400 ml))
- Distal Attachment for Direct Socket (Distal Attachment Pin and Distal Attachment Nut; D-040532)
- Taping Ring (D-040552)
- Lanyard Cord (L-392018)
- Relief Pads (M-100001)
- Icecast Case

INDICATIONS FOR USE

The Direct Socket Tool Kit is for direct lamination using the Direct Socket system, following trans-tibial amputation requiring a socket for a residual limb to be used with a prosthesis.

CONTRAINDICATIONS

None known.

INTENDED USE

Direct Socket Tool Kit is used for direct lamination using the Direct Socket system.

For lower limb application only.

Icecast Bladder is not to be used on amputees with the following conditions:

- Ultra-short trans-tibial with amputation performed directly under the tuberosity of the tibia
- Long residual limb exceeding bladder length
- Liner size above size 34

GENERAL SAFETY INSTRUCTIONS

PRECAUTIONS:

Icecast Bladder is sensitive to puncturing, especially when inflated and inverted. It can be easily damaged by improper handling of sharp objects such as scissors, knives and cured braids. Fingernails may also cause

damage during donning and doffing.

Ensure Icecast Bladder is not inverted when not in use to minimize risk of damage.

Excessive pressure can cause damage to the chambers; follow the pressure chart guidelines and avoid exceeding the stated values (**Fig. 3**).

Do not inflate Icecast Bladder to above 100 mmHg.

CLEANING

Use lukewarm soapy water or alcohol to clean the device on both the inside and outside.

INSTRUCTIONS FOR OPERATION

To produce consistent results with the Direct Socket Tool Kit, perform the following steps:

1. Connect Icecast Pump nozzle into Icecast Bladder valve. A proper connection is achieved when pump audibly clicks into place.
2. Inflate Icecast Bladder to a chamber pressure of 20 mmHg. If overinflated, release air by pressing black Air Release Button below pressure gauge.
3. Invert Icecast Bladder by rolling it back.
4. Attach inner end of Icecast Bladder to Distal Attachment (**Fig. 1**).
5. Keeping Icecast Pump aligned straight with long axis of the residual limb and knee, roll Icecast Bladder over the residual limb (**Fig. 2**).
6. If resistance occurs, slightly deflate the Icecast Bladder for easier donning while rolling on. Deflating slightly lessens excessive elongation and avoids excessive pressure build-up.
7. Inflate Icecast Bladder to correct pressure level (**Fig. 3**).

Direct Socket pressure casting	New amputee	Fleshy	Muscular	Bony
mmHg	40 – 50	60	70	80

8. Disconnect Icecast Pump from Icecast Bladder and wait 10 minutes for material to cure.
9. Deflate and remove Icecast Bladder from the residual limb (**Fig. 4**).

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

WARRANTY

Direct Socket Tool Kit comes with a 12-month warranty.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/ her clinical specialist. If any problems occur with the use of this product, immediately contact your medical professional.

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