ENGLISH

Note: Before sizing, please remove and discard the separator between the FRONT and Sternum collar pieces.

SIZING

Proper sizing is important for patient immobilization and comfort. Always maintain the patient's head in neutral alignment.

Neck Height

Place the chin-cup of the front to the patients chin and slide the lower front down to the patient's sternum (Figure 1). Lift the front away from the patient and push the two tabs to lock the size (Figure 2).

Warning: This collar is designed to be sized and locked into a single position only once; it is not intended to be resized.

Note: If you are familiar with the previous Philadelphia sizes you can use the size numbers marked on the front parts where 2 refers to 2 1/4", 3 refers to 3 1/4", etc. Simply slide the lower front down until the correct size number appears in the size window and lock.

Circumference

Identify your patient's profile type according to the chart in (Figure 3). Adjust the circumference by tearing away sections of the back piece foam to change to a smaller circumference (Figure 4).

Collar Sizing Guide (Neck height and circumferential adjustment):

The table below represents common sizing for the illustrated patient types (Figure 3). Actual sizing must be verified by the healthcare provider for each patient individually.

Patient Type	Neck Height	Circumference	Circumference Adjustment
Stout	2	Large	No tear off
Average male	3	Medium	Tear off 1 section on each side
Thin female	4	Small	Tear off 2 sections on each side

APPLICATION

- Proper application of the collar is as important as proper sizing. The two-piece design consists of a FRONT piece and a BACK piece which are packaged as a set. After adjusting the size of the collar, apply the back piece of the collar to the back of the patient's neck. Center the collar (Figure 5). The arrow on the back should point upward.
- Apply the FRONT piece of the collar with the chin secured in the recess. Center the collar to secure neutral alignment. The front piece OVERLAPS the back piece to ensure effective immobilization and comfort (Figure 6). The arrow on the front should point upward.
- With hook and loop fasteners, tighten the collar with a bilateral adjustment (Figure 7). This will secure the patient's cervical region in neutral alignment.

FEATURES

- Non-toxic Plastazote® Foam
- Adjustable to most patients for custom fit
- Radio Lucent
- Easily cleaned with soap and water
- Hook and loop fastener
- Latex-free

Note: Use as prescribed by a physician.

Caution: Federal law restricts this device to sale by, or on the order of, a physician.

Foam: Polyethylene
Frame and Rivets: Polypropylene
Straps: Nylon

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©Copyright Össur 2019