# Instructions for use

## Total Shock®

4400

4410

*Life Without Limitations®*

.....	3
ⒺⒺ Instructions for use.....	5
ⒺⒺ Gebrauchsanweisung.....	9
ⒺⒺ Notice d'utilisation.....	13
ⒺⒺ Instrucciones para el uso.....	17
ⒺⒺ Istruzioni per l'uso.....	21
ⒺⒺ Brugsanvisning.....	25
ⒺⒺ Bruksanvisning.....	28
ⒺⒺ Gebruiksaanwijzing.....	31
ⒺⒺ Instruções de Utilização.....	35
ⒺⒺ 取扱説明書.....	39

**EN - Caution:** Össur products and components are designed and tested to ISO 10328. Compatibility and compliance with this standard is achieved only when Össur products and components are used with other recommended Össur or authorized components. If un-usual movement or product wear is detected in a structural part of a prosthesis at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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




---

## Total Shock®

Two elastic control elements allow independent adjustment for compression and rotation. Compression is limited to 14mm of vertical displacement (**Figure A**). The Total Shock® allows for 13 degrees of in-rotation and 10 degrees of outrotation if installed according to instructions.

The Total Shock® is compact and can easily be incorporated into any endoskeletal system utilizing a 30mm tube. Trans-femoral (AK) and transtibial (BK) amputees can both benefit from this unique design.

## The Total Shock® model 4400 and the model 4410:

- Made of light weight aircraft alloy
- Model 4400: weighs 270 grams (9.5oz)
- Model 4410: weighs 300 grams (10.6oz)
- Weight limit: 100 kg (220lbs)

## Model 4400

Tighten 5mm bolt to 8 Nm (6 lb-ft) Dimensions:  
Length at Maximum Extension 139mm (5.47")  
Length at Maximum Compression 125mm (4.92")  
Maximum Travel 14mm (0.55")  
Diameter (at largest point) 44.0mm (1.73")

## Model 4410 (Figure B)

Tighten four 8mm set screws to 11 Nm (8 lb-ft)  
Dimensions:  
Length at Maximum Extension 146mm (5.74")  
Length at Maximum Compression 132mm (5.19")  
Maximum Travel 14mm (0.55")  
Diameter (at largest point) 44.2mm (1.74")  
Weight 300g (10.6oz)

## Vertical Compression Element Selection

These recommendations are meant as starting points only. You may have to move up or down the chart based on amputee's comfort, activity levels, or element compression.

COLOR	DUROMETER	PATIENT WEIGHT LBS & KG
Black	Regular	Approx. 88 - 143 lbs. (40 - 65kg)
Red	Firm	Approx. 143 - 198 lbs. (65 - 90kg)
Yellow	Extra Firm	Approx. 198 - 220 lbs. (90 - 100kg)
Orange	Extra, Extra Firm	Approx. 198-220 lbs. (90-100kg)

**Note:** models 4400 and 4410 are included with black, red and yellow elements. The Total Shock® is preinstalled for the right side with the pylon receiver or female pyramid installed distally. Inverted installation (pylon receiver or female proximally) will suit for the left side.

## Initial Fitting

When the amputee is standing, with equal weight on both limbs, the Total Shock® should compress a maximum of 1-3 mm. The Compression resistance can be adjusted by either changing the compression element or by turning the adjustment plug in the top of the Total Shock® with a coin.

## TORSIONAL FORCE REDUCTION

### (Torque Suppressor)

The Total Shock® also works as a torque absorber or suppressor. This function gives greater comfort to the amputee by reducing the shear

force produced between the proximal socket brim and the amputee's residual limb. This feature can be independently adjusted to accommodate the amputee's individual needs. The degree of torque absorption is dependent upon the level of amputation, design of socket, and corresponding socket fit.

### Rotational Control Element Selection (Figure C)

These recommendations are meant as starting points only. You may have to move up or down the chart based on amputee's comfort, activity levels or element compression.

COLOR	DUROMETER	PATIENT WEIGHT LBS & KG
Black	Regular	Approx. 88 - 143 lbs. (40 - 65kg)
Red	Firm	Approx. 143 - 198 lbs. (65 - 90kg)
Yellow	Extra Firm	Approx. 198 - 220 lbs. (90 - 100kg)
Orange	Extra, Extra Firm	Approx. 198-220 lbs. (90-100kg)

**Note:** models 4400 and 4410 are included with black, red and yellow elements.

### Initial Fitting

Have the amputee stand on a surface, which provides the highest possible coefficient of friction between the shoe and weight-bearing surface (**example:** carpet, asphalt, etc.) Now have them shift their weight onto the prosthesis and try to rotate the prosthesis as much as possible. The amount of rotation should be approximately 5 degrees for out rotation and more for in rotation. If it rotates too easily, you should change the rotational control element to one with a greater durometer or stiffness.

### Lubrication of Elastic Elements

The original elements have been lubricated (compression and rotation), however, if it is necessary to change an element, it must be lubricated to ensure proper function. Note: If the Total Shock® is installed with the retaining clamp positioned distally CentraLube™ oil should be used for element lubrication. If the Total Shock® is installed with the retaining clamp positioned proximally CentraLube™ grease should be used for element lubrication. A light but complete coverage of lubricant is required in either application.

Model 4401 CentraLube™ oil (Included).

Model 4401-G CentraLube™ grease (Not included) For inverted applications.

### HOW TO CHANGE THE URETHANE CONTROL ELEMENTS FOR COMPRESSION AND ROTATION (FIGURE D+E)

Be certain to keep all the parts absolutely clean when opening the Total Shock®.

Any kind of abrasive contamination on the inner surfaces will dramatically shorten the life of the product.

1. Remove the proximal Adjustment plugs by rotating it counterclockwise.
2. Remove the elastic compression element.
3. Remove the retaining ring from the end of the 8 mm axle with a small screwdriver. Always remove the proximal adjustment plug before removing the rotation control axle.
4. Carefully remove the 8mm axle containing the rotational control element by pushing it through the large opening. **(Figure F)**
5. Slide off the rotational control element. Apply a light, complete coating of CentraLube™ oil or grease (depending on shock orientation) to the new element and replace.

6. After applying a light, complete coating of CentraLube™ oil or grease (depending on shock orientation) to both the proximal and distal ends of the compression element, place the element into the proximal opening of the Total Shock®. Reinstall the proximal adjustment plug by rotating clockwise.
7. Installing the perforated rotational elements: The holes shall face posterior on the lateral side and anterior on the medial side to allow for more stability. **(Figure G)**
8. Turn the adjustment plug **CLOCKWISE** for stiffer compression. **COUNTERCLOCKWISE** for a softer compression.

**Note:** The minimal setting of the proximal adjustment plug should be flush with the top edge of the Total Shock®.

### **CAUTION!!**

The Total Shock® is designed to be used in conjunction with an aluminum or fiber composite tube with insert.

**Please be advised:** If you are using a Total Knee model 1900, 2000 or 2100 with the Total Shock® model 4400 or 4410 and you install the extension assist model 2070 or 2070S you should use the A-433120 (100kg), A-435120 (166kg) Long Distal pyramid Adaptor connected to the A-342100 (100kg), A-345100 (166kg) Tube Clamp Adaptor. Either one of the A-433110 (100kg), A-435110 (166kg) or A-433120 (100kg), A-435120 (166kg) can be used in the distal end of the Total Shock® model 4400.

The clamp must be positioned posteriorly (never medially or laterally). **(Figure H)**

Orange Vertical Compression Element is optional in models 4400 and 4410 when the Yellow Compression Element does not provide sufficient resistance. This may be appropriate for your high impact use, i.e. recreational activities.

### **Maintenance and Inspection**

The Total Shock® is designed and manufactured to provide long and trouble free service. To ensure proper function it is recommended that the Total Shock® be inspected for unusual signs of wear, cleaned and lubricated every six months.

#### **Six Months Inspection**

Visual inspection checklist:

1. No cracks in housing.
2. Vertical compression and rotational control elements are not worn or cracked.
3. Large and small O-rings on axle are in place and show no signs of damage or wear.
4. Rotational control element is in position with its widest side horizontal.
5. Retaining ring is in place.

#### **Cleaning**

In a clean environment completely disassemble the Total Shock®. Wipe all surfaces of the Total Shock® with a soft cloth moistened with a small amount of kerosene or acetone. **Note:** Do not use any stronger solvents.

#### **Lubrication**

**Note:** If the Total Shock® is installed with the retaining clamp positioned distally CentraLube™ oil should be used at all lubrication points. If the Total Shock® is installed with the retaining clamp positioned proximally CentraLube™ grease should be used at all lubrication points. Once the components of the Total Shock® are clean, apply a light, complete coating of the prescribed CentraLube™ lubricant to the two internal white bearing surfaces. Reassemble the piston and housing. Apply a

thin, complete coating of the prescribed lubricant to all exterior surfaces of the rotational element and install. Apply a light, complete coating of the prescribed lubricant to both the proximal and distal ends of the linear element and install. Wipe down all external surfaces and install in prosthesis.

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### **Össur Iberia S.A.**