Instructions for use

Unloader<sup>®</sup> Spirit

*Life Without Limitations<sup>®</sup>*

.....	3
ⓔⓃ Instructions for use.....	5
ⓓⓔ Gebrauchsanweisung.....	8
ⓕⓗ Guide de fabrication.....	11
ⓔⓢ Instrucciones para el uso.....	14
ⓔⓣ Istruzioni per l'uso.....	17
Ⓝⓞ Bruksanvisning.....	20
ⓓⓐ Brugsanvisning.....	23
ⓈⓋ Bruksanvisning.....	26
ⓕⓙ Käyttöohjeet.....	29
Ⓝⓕ Gebruiksaanwijzing.....	32
ⓔⓟ Instruções de Utilização.....	35

**EN - Caution:** This product has been designed and tested based on single patient usage and is not recommended for multiple patient use. If any problems occur with the use of this product, immediately contact your medical professional.

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
## INDICATIONS FOR USE

- Mild to Severe Unicompartmental Osteoarthritis
- Unicompartmental knee conditions that require load reduction to affected compartment such as articular cartilage defect repair, meniscal cartilage repair, avascular necrosis, and tibia plateau fracture
- Normal activities of daily living
- Patients under 250 lbs (115Kg)

## CAUTION

For patients with peripheral vascular disease, neuropathy, and sensitive skin, close physician supervision is recommended.

The Unloader®Spirit knee brace is modifiable to provide a good fit to the leg anatomy while providing corrective forces. Modifications to extension stops, shell alignment, varus or valgus force, strap length, pad size, and shell contour should only be made by the health care professional.

### Adjusting hinge contour and corrective forces (Picture 1):

1. Before putting the brace on the leg, loosen the two tamper-proof hinge adjustment screws using the supplied Allen wrench.
2. Seat the patient on the edge of a chair, leg in full extension
3. Loosen the tibia strap
4. Place brace on leg with center of hinge in line with center of patella. The hinge should be placed midline in the anterior / posterior plane
5. Secure midstrap. Use spacer pads to quantify hinge angulation (each pad = 3 degrees). (Picture 2).
6. Insert spacer pads between the patient's leg and the condyle pad of the hinge
7. Align shells and hinge with leg anatomy
8. With slight pressure against the hinge, tighten the bottom hinge adjustment screw followed by the top adjustment screw
9. Remove the brace, tighten the adjustment screws
10. To ensure adjustment is secure, apply varus or valgus strain to the mechanism. If the mechanism clicks, go to step 1

### Shell adjustment:

1. Remove liners from the brace
2. Use heat gun to warm specific area of shell
3. Apply pressure to adjust the shell shape

### Strap length adjustment (Picture 3):

1. Determine strap length to be removed
2. Remove strap from brace
3. Locate Velcro loop end and remove Velcro hook piece
4. Cut along sewn cut lines. Each line is 1 inch (2.5 cm)
5. Re-apply strap

### Strap pad adjustment (Picture 4):

When strap pad interferes with D-ring or ability to tighten strap:

1. Remove pad from strap
2. Cut along indentations
3. Place back on strap

### **Applying the brace for use:**

*Note: Instructions show medial hinged brace*

1. Sit patient on edge of chair, leg fully extended
2. Loosen the tibia strap (**Picture 7**)
3. Position brace on leg with center of hinge in line with top of kneecap  
Attach the brace straps in the following sequence.
4. Mid Strap (**Picture 5**)
  - a. Wrap strap around back of leg with strap pad on top of calf muscle
  - b. Thread through strap ring and secure
  - c. Straps should be tight, without restricting circulation or causing cramping
5. Calf Strap (**Picture 6**)
  - a. Wrap strap around calf
  - b. Thread through the strap ring and secure
  - c. Do not over tighten
6. Tibia Strap (**Picture 7**)
  - a. Tighten Tibia strap to be snug, yet comfortable
  - b. Do not pull the shell out of position
  - c. The strap does not need adjustment each time the brace is applied
7. Thigh Strap (**Picture 8**)
  - a. Wrap strap around back of thigh - angle slightly down
  - b. Thread through the strap ring and secure
  - c. During use, the thigh strap may need to be re-tightened
8. The Dynamic Force Strap™ (DFS) (**Picture 9**)
  - a. Wrap strap around the back of leg
  - b. Thread through D-ring ring on the thigh shell
  - c. Bend braced knee to an 80-degree angle
  - d. Hold thigh shell or hinge to prevent brace from moving
  - e. Tighten DFS. The DFS should be snug
  - f. Do not over tighten DFS, as this pulls the brace out of position

### **Range of Motion Adjustment (ROM)(Picture 10):**

The brace comes with a 0° extension stop setting. The extension stops can be adjusted using the 5°, 10°, or 15° stop. *Note: Optimal corrective forces are applied at full extension.*

1. Remove the condyle pads from the hinge
2. Remove retaining screw
3. Slide the 0° stop out of the hinge
4. Insert desired extension stop
5. Secure with the retaining screw. **DO NOT OVER TIGHTEN!**

### **Final Inspection: Right leg braces (Pictures 11 & 12, ):**

11. Medial Front
12. Medial Back

### **BRACE MAINTENANCE**

The Unloader Spirit brace is designed for minimum maintenance and care. Follow these easy maintenance steps to keep the brace in top condition.

#### *Cleaning*

Liners, Straps and Strap Pads:

- The liners mid strap, calf strap, thigh strap and all strap pads may be removed for cleaning
- **DO NOT REMOVE THE DFS STRAP**
- **DO NOT MACHINE WASH. DO NOT PLACE IN DRYER**
- Hand-wash using mild non-detergent soap

- Rinse thoroughly
- Air-dry

#### Hinge:

- For safety reasons, the hinge should only be disassembled by a health care professional
- Check regularly for foreign materials (ie. dirt or grass)
- Clean dirt from hinge using tap water or compressed air

If the brace is used in salt water or chlorinated water, rinse well and air dry brace. A regular service routine will ensure the hinge, liners, straps and pads are in optimal working order.

#### **BRACE ACCESSORIES**

Ossur offers accessories that meet the specific needs of the Unloader® Spirit brace. These accessories can be purchased through the brace provider or Ossur Customer Service Departments.

- Extra Set of Liners
- Replacement Strap Set
- Undersleeves for use with sensitive / thin skin
- Patella Control Sleeve
- Sports Cover
- All Purpose Bags

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Unloader® Spirit is covered by one or more of the following patents. Other U.S. and foreign patents pending. ® indicates trademark registration in U.S.A. and selected countries only.