




Instructions for Use

---

CHEETAH  
(FLEX-SPRINT III)

|    |                           |    |
|----|---------------------------|----|
|    |                           | 3  |
| EN | Instructions for Use      | 5  |
| DE | Gebrauchsanweisung        | 8  |
| FR | Notice d'utilisation      | 11 |
| ES | Instrucciones para el uso | 14 |
| IT | Istruzioni per l'uso      | 17 |
| DA | Brugsanvisning            | 20 |
| SV | Bruksanvisning            | 23 |
| NL | Gebruiksaanwijzing        | 25 |
| PT | Instruções de Utilização  | 28 |
| JA | 取扱説明書                     | 31 |

# ENGLISH

---

## CHEETAH SYSTEM

- Custom foot for track & field sports
- Preferred transtibial sprinting foot
- Clearance for hardware  
Minimum: 250mm (10") for Pylon Adapter (D)  
Maximum: 32mm (1 1/4") for lamination adapter (C)
- Maximum weight 147 kg (325lbs) for lamination adapter
- One adult foot size

## CONNECTOR OPTIONS (A)

Pylon Connector FSX50003 (1-6), FSX50004 (7-9)

- Used for clearance greater than 250mm (10")
- Adapter is laminated direct to the back of the socket protruding from the apex of the distal socket to allow attachment of the foot.

Bolts go through adapter first then through foot module. Apply Loctite 410 to all threads and torque to following specifications:

| Category | Bolt Size | Torque            |
|----------|-----------|-------------------|
| M 1-6    | 5/16"     | 33 Nm (24 ft-lbs) |
| L 7-9    | 3/8"      | 45 Nm (33 ft-lbs) |

## LAMINATION CONNECTOR (B)

FSX50001 (1-6), FSX50002 (7-9)

- Used for clearance less than 250mm (10")
- Adapter is laminated direct to the back of the socket

Bolts go through adapter first then through foot module. Apply Loctite 410 to all threads and torque to following specifications:

| Category | Bolt Size | Torque              |
|----------|-----------|---------------------|
| M 1-6    | 5/16"     | 33 Nm ( 24 ft-lbs ) |
| L 7-9    | 3/8"      | 45 Nm ( 33 ft-lbs ) |

**THE CHEETAH™ MAY ALSO BE CONNECTED BY A DIRECT LAMINATION.**

## FRICION PADS

Always use between T and foot module. If alignment wedges are used, place friction pads between each layer.

## BENCH ALIGNMENT

### Sagittal plane (C+D)

- Load line bisects socket and falls through the ground contact point or approximately 5cm (2") posterior.
- Determine the appropriate socket flexion according to amputee requirements.

### Coronal plane

- Introduce appropriate socket abduction/adduction (using a jig)
- 5°-7° additional external rotation
- Achieve proper height approximately 25mm (1") longer than the

sound side with shoe to compensate for vertical compression and plantar flexion

**Note:** The recommendations for bench alignment are used as a starting point. Definitive alignment will be achieved during dynamic assessment.

### **CUTTING FOOT MODULE**

- Before cutting check for enough clearance.
- Measure for length and add 50mm (2") to compensate for measurement error. The final height is approximately 25mm (1") taller than the sound side.
- Clamp foot module in vice. Use a sharp fine tooth hacksaw blade.
- Sand away any sharp edges.
- Secure foot to using fitting clamps or synthetic plaster bandages.
- Cut again if necessary.
- Lowest hole position is in 30 cm (12") when pylon is vertical.

### **DRILLING FOOT MODULE ATTACHMENT HOLES**

- Back-up foot module when drilling using wood block or discarded carbon fiber section and clamp securely.
- Use drill guide template to locate holes.
- Use sharp drill bits. Suggestion: Bullet nose or solid carbide tipped drill bits.
- Drill size:
  - Cat 1-6 8,4mm (21/64")
  - Cat 7-9 11mm (7/16")
- Drill speed 1200-1800 RPM.
- Feed drill slowly with light intermittent pressure - do not force.
- Foot is now ready for assembly.

### **DYNAMIC ALIGNMENT**

Physical ability and length of residual limb will result in different alignment. Best performance will be achieved with best achievable symmetry of gait. Each athlete will have her/ his own optimal alignment, which relates to the individual's running style. Optimal alignment is critical in order to reduce energy expenditure and eliminate excessive muscle and ligament strain while running.

- Check for best symmetry of gait
- Stride length can be affected by A-P position socket to foot
- Vertical deflection is in the average between 50mm (2") to 75mm (3")

### **TROUBLE SHOOTING**

Exclude reasons related to socket fit. Try to achieve symmetry in running as close as possible. The symmetry in running should be assessed at maximum speed.

### **FURTHER RECOMMENDATIONS:**

For track & field events shoes usually feature spikes:

The sole of the spare shoe is then attached to the prosthetic foot in the following way:

- Cut off the sole area where the spikes are attached (1).
- Remove the insole
- Bond leather to the back side of the prepared part.
- Sand the surfaces (2).
- Bond assembly to the ground contact area (3).

**WATERPROOF**

A waterproof device allows use in a wet and/or humid environment and temporary (<30 mins.) submersion of up to 1m in fresh water. Exposure to fresh water under these circumstances shall have no harmful effect.

Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

The components used with the foot need to be waterproof as well.

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**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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