ENGLISH

CONTENT

- Lock Body with guide (w/o Dummy Screw) (**fig. 3**).
- Dummy Screw (**fig. 4**).
- Clutch Mechanism (**fig. 7**) and push button (**fig. 8**).
- T-wrench (**fig. 10**).
- Attachment Pin clutch, standard (**fig. 11**).
- Delrin Guide (**fig. 6**).
- Silicone Sealing (**fig. 5**).
- Pin Tunnel w/ o-ring (**fig. 2**).
- Pyramid Bolt long w/ o-ring (**fig. 1**).

Specially designed Pyramids for single lamination are available:

L-273000 Icelock Stainless Pyramid 273 (rated up to 100kg/220lbs).

L-272000 Icelock Titanium Pyramid 272 (rated up to 166kg/365lbs).

LAMINATION INSTRUCTIONS

- Flatten distal end of Plaster Cast (9) to center Delrin Guide (6). Attach with nails provided or double-sided tape maintaining original pin alignment.
- Apply inner PVA bag and tie it around neck (G) on Delrin Guide (6).
- Trim away excess PVA and centrally locate Silicone Seal (5).
- Place Lock Body (3) onto Delrin Guide (6), ensuring its rotation allows for correct positioning of Clutch Mechanism (7).
- Secure Lock Body (3) and Pin Tunnel (2) with assembly bolt (1) (making sure that o-ring is between bolt / Pin Tunnel and Pin Tunnel / Lock Body) using 5mm hexagon key.
- NOTE: Be sure to tighten so that resin will not leak in between components.
- Place some clay in hole on bolt (A) and groove (B) on Dummy Screw (4).
- Tie half of lamination material into groove on Lock Body (C) and remainder over end of Lock Body (D) leaving Pin Tunnel (2) and bolt head (E) exposed.
- Apply outer PVA bag and finish lamination.
- When curing has finished, remove bolt (1), expose but do not remove Dummy Screw (4.)
- Prior to second lamination ensure that Pin Tunnel (2) is sealed.
- After completion of second lamination, expose and remove Dummy Screw (4).
- Insert and tighten Clutch Mechanism (7) to 15Nm (11ft.lb) - VERY IMPORTANT.

NOTE: Push Button (8) may be shortened by unscrewing from lock body and cutting it with hacksaw.

THERMOPLASTIC INSTRUCTIONS

- Flatten distal end of Plaster Cast (9) to center Delrin Guide (6). Attach with nails provided or double-sided tape maintaining original pin alignment.
- Place Lock Body (3) onto Delrin Guide (6), ensuring its rotation allows for correct positioning of Clutch Mechanism (7).
- Secure Lock Body (3) and Pin Tunnel (2) with assembly bolt (1) using a 5mm hexagon key.
- Vacuum form thermoplastic, ensuring that plastic conforms well to Lock Body (3) and Pin Tunnel (2).
- Once plastic has cooled, trim around Pin Tunnel (2) and Dummy Screw (4). IMPORTANT: Before making outer socket, always ensure that assembly bolt (1) has been removed and Pin Tunnel (2) resealed.
- After manufacture of outer socket, expose and remove Dummy Screw (4.)
- Insert and tighten the Clutch Mechanism (7) to 15Nm (11ft.lb) - VERY IMPORTANT.

NOTE: Push Button (8) may be shortened by unscrewing from lock body and cutting it with hacksaw.

Attachment Pin: Apply Loctite 410 or 411 onto threads (F) on Attachment Pin (11) and screw it into Locking Liner. Torque Attachment Pin to 4Nm (3ft.lb).

Also available:

L-292003 Attachment Pin clutch, short (7 notches).

L-292005 Attachment Pin clutch, long (17 notches).

NOTE: Use attachment pins for Icelock 200 series only.

T-wrench (10) supplied is specially designed to wind Push Button (8) on Clutch Mechanism (7) when donning socket. Use of any other tool may result in over tightening.

- Limited one-year warranty.
- Rated up to high activity level (K3).

A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

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EN – Caution: Ossur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Ossur products are used with other recommended Ossur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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Ossur Deutschland GmbH