




Instructions for Use

AIRFORM[®] UNIVERSAL
INFLATABLE STIRRUP
ANKLE BRACE

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	5
FR	Notice d'utilisation	6
ES	Instrucciones para el uso	7
IT	Istruzioni per l'uso	8
NO	Bruksanvisning	9
DA	Brugsanvisning	10
SV	Bruksanvisning	11
FI	Käyttöohjeet	12
NL	Gebruiksaanwijzing	13
PT	Instruções de Utilização	14
PL	Instrukcja użytkowania	15

ENGLISH

USES

Mild to moderate sprains, strains, post cast support and injury prevention.

CAUTION

Patients should be instructed to inspect their ankle, foot, and product before and after each use. Always wear an absorbent sports sock.

PHOTO CAPTIONS

- A. Adjustable contact closure straps
- B. Side stirrups
- C. Adjustable heel strap
- D. Air tube
- E. Air valve
- F. Air bladder

FITTING INSTRUCTIONS

1. Step into the brace to check fit. Center the side stirrups along lower leg and ankle.
2. Adjust heel width, if necessary, by lifting the lower air bladder or foam pad and adjusting the contact closure. The heel strap should be adjusted equally on both sides. Position the bottom of the stirrups to be at or slightly above the bottom of the foot.
3. Wrap lower strap snugly around brace.
4. Wrap top strap and readjust the lower strap as necessary to achieve a snug fit.
5. Apply shoe and tighten laces. Readjust the brace as necessary for maximum comfort and support.
6. If necessary, inflate air bladders for additional compression and support. Inflate by peeling back the top portion of the bladder to locate the air valve. Insert the provided air tube through the air valve and inflate. Do not overinflate. To remove air, place tube into air valve and press on the air bladders until desired inflation is achieved.

WASHING INSTRUCTIONS

This support may be cleaned by wiping with a damp cloth and mild detergent. Do not use bleach or strong detergents. Do not wash in washing machine or place in dryer. Avoid direct contact for prolonged time in sun or direct heat.

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