<ul style="list-style-type: none">Material: Titanium. Weight Limit: 166kg. Weight: 112g.		
	Torque	Loctite
Mounting Screws	10 Nm	243

For secure fastening use the loctite recommended and apply rated torque twice.
May be used for Transtibial and Transfemoral laminated sockets.

INSTRUCTIONS FOR INTEGRATION INTO LAMINATED SOCKETS.

Socket must be laminated in two layers. The 4-Hole Socket Adapter should be attached to the inner layer with putty. All surfaces should be roughened prior to second lamination.
1. Isolate cast with PVA bag.
2. Apply appropriate lay up material for first lamination.
3. Apply outer PVA bag and laminate.
4. Once cured, remove PVA bag and roughen outer surface of lamination.
5. Attach 4-Hole socket adapter in correct alignment and rotation with putty.
6. Position Lamination Cap on Adapter and secure using a 4mm Hex key.
7. Seal putty and Lamination Cap with sealing resin.
8. Apply Lay up material and combinations of materialappropriate to the individuals weight and impact/activity level.

IMPORTANT: Tie 50mm carbon reinforcements into both grooves, between lamination cap and adapter, and groove in adapter

The Lamination Cap should be removed following lamination by unscrewing the bolts using a 4mm Hex key.

- 2 year limited warranty.**

Titanium tested to ISO 10328, load level A100+66%, rated to 366lbs/166kg for normal use.

A waterproof device allows use in a wet and/or humid environment and temporary (<30 mins.) submersion of up to 1m in fresh water. Exposure to fresh water under these circumstances shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

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