



MIAMI JR® CERVICAL COLLAR

- Innovative collar designed exclusively for children 12 and under
- Four true pediatric sizes ensure optimal airway and the upper cervical spine immobilization
- Sized according to the Broselow-Luten Color-Coding System
- Sorbatex™ pads make the collar “skin friendly” and protect skin integrity during extended wear
- The Össur pediatric system also includes the NecLoc® Kids Extrication Collar, Occian™ AirWay PAD™ for appropriate positioning and the Lil'Angel MR-Conditional (Safe at 3-Tesla) Halo System

| COLLAR | COLLAR SET | REPLACEMENT PADS | SIZE |
|--------|------------|------------------|--------------|
| MJ-P0 | MJR-PO | MJP-PO | 0-6 mths |
| MJ-P1 | MJR-P1 | MJP-P1 | 6 mths-2 yrs |
| MJ-P2 | MJR-P2 | MJP-P2 | 2-6 yrs |
| MJ-P3 | MJR-P3 | MJP-P3 | 6-12 yrs |

Occian AirWay Pad (part #AP-T10) available.

Replacement pads and collars may be ordered separately.




PAPOOSE®

- Securely and gently immobilizes the head and spine in the proper anatomical airway and spinal alignment
- Occipital shape protects against positional plagiocephaly
- Allows for continued intimacy between caregiver and child
- True cervical spine immobilizer for infants up to 3 months

| PART# | DESCRIPTION | SIZE |
|---------|--------------------------------|------------------------|
| PO-100 | Papoose without P0 Front | Infants up to 3 months |
| PA-100 | Papoose with P0 Front | Infants up to 3 months |
| PAR-100 | Papoose with P0 Front and pads | Infants up to 3 months |

Additional replacement pads available



MIAMI JTO THORACIC EXTENSION

- Provides more control than a collar alone can provide
- Innovative device designed to treat cervical and high-thoracic injuries
- Patented design with floating back permits maximum immobilization without the risk of compensatory motion
- Snap on and off without any tools
- Sorbatex padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance

| PART# | DESCRIPTION |
|---------|---------------------------|
| MT-900 | JTO |
| MTR-900 | JTO with replacement pads |
| MTP-950 | Extra replacement pads |