



## PHILADELPHIA® TRACHEOTOMY COLLAR

- Anterior and posterior rigid plastic reinforcement supports limit movement
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and aqua-therapy
- Easily customized for individual requirements
- Not made with natural rubber latex
- Non-toxic and hypoallergenic Plastazote® foam reduces skin irritation

PART#	HEIGHT	SIZE	CIRC.
PHP-TI	1.25"	Infant	6-8"
PHP-TP	1.75"	Pediatric	8-11"
PHP-TXXS	see footnote	XSmall	up to 10"
PHP-TYS	see footnote	Small	10-13"
PHP-TYM	see footnote	Medium	13-16"
PHP-TYL	see footnote	Large	16-19"
PHP-TYXL	see footnote	XLarge	19" up

X = 2 (2.25" height)

Y = 2 (2.25" height); 3 (3.25" height); 4 (4.25" height); or 5 (5.25" height)

### Philly Liners

PART#	DESCRIPTION
PHP-LCS	All XSmall/Small sizes
PHP-LCM	All Medium Sizes
PHP-LCL	All Large Sizes
PHP-LCXL	All XLarge Sizes