

AFO LIGHT

AFO Light provides support only where it is needed, eliminating unnecessary bulk and providing a cosmetically appealing look. Using the best of Flex-Foot® technology and the energy storing properties of carbon fiber, AFO Light offers a strong and lightweight solution for people with mild to moderate drop foot.

- Energy return and propulsion assistance of the limb at toe off
- Supports the foot for ground clearance during swing phase
- Carbon fiber design provides excellent strength-to-weight ratio and dynamic motion
- Gradual plantar flexion at initial contact
- Stiffness gradient at the heel and the toe correspond to ground reaction forces as they are applied to the foot
- Easy to don and doff

INDICATIONS

- Mild to moderate drop foot
- Contra-indicated for severe ankle-foot deformities, severe spasticity, fluctuation edema, M-L instability, or extreme activity
- Suggested weight range: up to 220 lbs (100kg)

	SMALL	MEDIUM	LARGE	XLARGE
Part number	O-1210X2	O-1210X4	O-1210X6	O-1210X8
Shoe size				
• Male	Up to 7.5	8-10.5	11-14	12-15
• Female	7-9	9.5-12		
Foot size	8.75-9.5" (22-24cm)	9.5-10.5" (24-27cm)	10.5-11.5" (27-29cm)	11.5-12.5" (29-32cm)
Foot plate	8.75" (22cm)	9.25" (23.5cm)	10.25" (26cm)	10.75" (27.5 cm)
Calf height¹	12.5-14.25" (32-36cm)	13.5-15" (35-38cm)	14.25-15.75" (36-40cm)	15-16.5" (38-42cm)
Spring length²	13.5" (34.5cm)	14.25" (36cm)	15" (38cm)	15.75" (40cm)
Heel height	0.6" (1.5cm)	0.6" (1.5cm)	0.6" (1.5cm)	0.6" (1.5cm)
Weight	2.8oz	3.15oz	3.5oz	3.85oz

X = 1 (Right Leg) or 2 (Left Leg)

¹ The height at the thickest part of the calf

² The length of the posterior spring

PART# ACCESSORY

O-122101 AFO Light Accessory Kit*

*AFO Light Accessory Kit includes strut padding, strut hook-and-loop, shrink tube, calf piece, calf piece padding, calf piece light hook-and-loop, light strap, light elastic band and silicone padding kit.

