



## REHAB

The Contour is designed to fit specifically the right or left leg, and can be ordered in either a long or short model. The Universal model is standard to fit both left and right. Both are available with breathable “Air Light” padding or tissue controlling “Full Foam”.

- Patient comfort with the cool and breathable Air Light design
- Warmth and maximum tissue control via Full Foam model.
- Yellow strap tabs allow quick brace placement and adjustment
- Push-button system for easy length adjustment
- ROM control in 15° increments: Flexion 0°-105°, Extension 15°-90°
- Low profile: fits easily under clothing
- Drop-lock option is available for some models

## INDICATIONS

- For post-operative and post-injury rehabilitation
- Suggested weight range: up to 250 lbs

### Rehab Contour

RIGHT#	LEFT#	FOAM	SIZE	DROP LOCK
B-220500000	B-220600000	Air Light	Regular	No
B-226500000	B-226600000	Air Light	Short	No
B-224500000	B-224600000	Air Light	Petite	No
B-221500000	B-221600000	Full Foam	Regular	No
B-233500000	B-233600000	Air Light	Regular	Yes
B-234500000	B-234600000	Full Foam	Regular	Yes

### Rehab Universal

PART#	LEG	FOAM	SIZE	DROP LOCK
B-222000000	Universal	Air Light	Regular	No
B-223000000	Universal	Full Foam	Regular	No
B-236000000	Universal	Air Light	Regular	Yes
B-235000000	Universal	Full Foam	Regular	Yes

### Replacement Softgoods

PART#	FOAM
B-705000502	Air Light
B-705000503	Full Foam