



## SHOULDER ABDUCTION SLING

- Provides shoulder abduction at 10° with support
- Pillow portion rests on waist to support weight of arm
- Snap-on sling design for easy application
- Easy strap length adjustability
- Breathable cotton provides a comfortable fit
- Fits left and right

PART#	SIZE	FITS ARM LENGTH*	U/M
204053	Small	up to 12.5"	1 ea
204055	Medium	12.5-15.5"	1 ea
204057	Large	15.5" and up	1 ea

\* Measure elbow crease to fingertip



## BUCKLE CLOSURE ARM SLING

- Contact closure shoulder strap with swathe

PART#	SIZE	LENGTH	DEPTH	STRAP	U/M
204020	Child	10.5"	5"	24"	6 pk
204022	XSmall	12.25"	4"	32"	6 pk
204023	Small	14.5"	5.25"	35"	6 pk
204025	Medium	17.5"	6.5"	39.5"	6 pk
204027	Large	19"	7.5"	43.5"	6 pk

## STRAP STYLE SHOULDER IMMOBILIZER

(photo not shown)

- Contact closure shoulder strap with clip buckle swathe

PART#	SIZE	LENGTH	DEPTH	STRAP	U/M
204013	Small	16"	6"	35"	1 ea
204015	Medium	17"	7"	39"	1 ea
204017	Large	18"	9"	42"	1 ea
204018	XLarge	19.5"	10"	46"	1 ea