


# ICEROSS DERMO SEAL-IN® LINER



Made with the softest silicone, Iceross Dermo cushions the limb while actively caring for the skin. Ideal for vascular patients and those with sensitive skin.

Iceross Dermo is a premium suspension liner that provides superb softness, gentle skin contact and unique Active Skin Care ingredients. The combination of DermoGel® silicone with an ultra-strong and elastic Supplex® fabric cover offers excellent durability and an intimate fit.

Össur recommends that Iceross Dermo Seal-In Liner is used in conjunction with the Icelock® 500 Series Expulsion Valve.



USER INFORMATION	
Amputation Level:	Transtibial
Impact Level:	Low to High
LINER INFORMATION	
Size Standard:	18, 20, 22, 23, 5, 25, 26, 5, 28, 30, 32, 34, 36.
Profile:	3mm
Matrix:	10cm

## SUSPENSION METHOD



SEAL-IN

## FEATURES



MATRIX



FABRIC COVER



ACTIVE SKIN CARE



SILKEN INNER SURFACE

## PART NUMBER INFORMATION - ICEROSS DERMO SEAL-IN LINER

PART#	PROFILE	SUSPENSION METHOD	MATRIX
I-4613XX	3mm	Seal-In	10cm