


ICEROSS SEAL-IN® X5 TRANSTIBIAL LINER WITH WAVE



The Iceross Seal-In X5 with Wave incorporates a series of five integrated seals that conform to the shape of the residual limb and the internal socket wall, providing an airtight seal.

The Iceross Seal-In X5 with Wave offers uncompromised security and freedom of movement.

Össur recommends that Iceross Seal-In X5 Transtibial Liner with Wave is used in conjunction with the Icelock® 500 Series Expulsion Valve.



USER INFORMATION	
Amputation Level:	Transtibial
Impact Level:	All
LINER INFORMATION	
Size Standard:	18,20,22,23.5,25,26.5,28,30,32,34,36.
Profile:	3mm
Matrix:	10-13cm

SUSPENSION METHOD




SEAL-IN X5

FEATURES



MINIMUM RECOMMENDED LENGTH OF RESIDUAL LIMB

LINER SIZE	MATRIX LENGTH
18	10cm
20-25	11cm
26.5-30	12cm
32-36	13cm



Matrix size varies depending on the liner size, the smaller the liner the shorter the matrix.

PART NUMBER INFORMATION - ICEROSS SEAL-IN X5 TRANSTIBIAL LINER WITH WAVE

PART#	PROFILE	SUSPENSION METHOD	MATRIX
I-3664XX	Wave 3mm	Seal-In X5	10 -13cm


ICEROSS SEAL-IN® X5 TRANSTIBIAL LINER



The Iceross Seal-In X5 incorporates a series of five integrated seals that conform to the shape of the residual limb and the internal socket wall, providing an airtight seal.

Ossur recommends that Iceross Seal-In X5 Transtibial Liner is used in conjunction with the Icelock® 500 Series Expulsion Valve.

USER INFORMATION	
Amputation Level:	Transtibial
Impact Level:	All
LINER INFORMATION	
Size Standard:	18,20,22,23.5,25,26.5,28,30,32,34,36.
Profile:	3mm
Matrix:	10-13cm



ICEROSS®
Transtibial – Seal-In®

SUSPENSION METHOD



SEAL-IN X5

FEATURES



MATRIX



FABRIC COVER




ACTIVE SKIN CARE



SILKEN INNER SURFACE

MINIMUM RECOMMENDED LENGTH OF RESIDUAL LIMB

LINER SIZE	MATRIX LENGTH
18	10cm
20-25	11cm
26.5-30	12cm
32-36	13cm



Matrix size varies depending on the liner size, the smaller the liner the shorter the matrix.

PART NUMBER INFORMATION - ICEROSS SEAL-IN X5 TRANSTIBIAL LINER

PART#	PROFILE	SUSPENSION METHOD	MATRIX
I-3663XX	3mm	Seal-In X5	10 -13cm