

# ICEROSS® SPORT LOCKING LINER




The ultimate suspension liner for active people.

For active amputees that experience greater stress on the residual limb and require additional stability and intimacy of fit. Iceross Sport features a unique double-layer silicone design and an Advanced Rotation Control System that together provide effective shock-absorption and rotational stability.

Unlike other suspension liners that are designed mainly for lower activity levels, Iceross Sport is designed specifically to give active amputees the stability and support needed to lead an active lifestyle.

Össur recommends that Iceross Sport Liner is used in conjunction with the Icelock® 600 Series.



USER INFORMATION	
Amputation Level:	Transtibial
Impact Level:	High to Extreme
LINER INFORMATION	
Size Standard:	16, 18, 20, 22, 23.5, 25, 26.5, 28, 30, 32, 34.
Profile:	3mm
Matrix:	10cm or Custom 2-14cm*

## SUSPENSION METHOD



LOCKING

## FEATURES



MATRIX



FABRIC COVER



ACTIVE SKIN CARE



DOUBLE DURO



MENTHOL

## PART NUMBER INFORMATION - ICEROSS SPORT LINER

PART#	PROFILE	SUSPENSION METHOD	MATRIX	ITEMS
I-1033XX	3mm	Locking	10cm	N/A
I-1233XX	3mm	Locking	Custom 2-14cm*	N/A
I-1033XXD	3mm	Locking	10cm	Double Pack

\*For a custom matrix, please specify the required matrix length on your order. Please note that deviation in the stabilizing matrix length can be +/- 1.5cm.