

ICEROSS DERMO® CUSHION LINER WITH WAVE




Iceross Dermo Cushion Liner with Wave is made with the softest silicone that actively cares for the skin.

This liner provides superb softness, gentle skin contact and Össur's unique Active Skin Care ingredients. The combination of DermoGel® silicone with an ultra-strong and elastic Supplex® fabric cover offers excellent durability and an intimate fit.

The Wave feature makes the liner easier to flex at the knee, allowing it to adapt more to the user's movements. Less resistance means that less energy is required to flex the knee.

Össur recommends that Iceross Dermo Cushion Liner with Wave is used in conjunction with Iceross Sleeve and Icelock® 500 Series.



USER INFORMATION	
Amputation Level:	Transtibial
Impact Level:	Low to High
LINER INFORMATION	
Size Standard:	16, 18, 20, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 40.
Profile:	3mm

SUSPENSION METHOD



CUSHION

FEATURES



FABRIC COVER



ACTIVE SKIN CARE



WAVE



SILKEN INNER SURFACE

PART NUMBER INFORMATION - ICEROSS DERMO CUSHION LINER WITH WAVE

PART#	PROFILE	SUSPENSION METHOD
i-cw43XX	3mm	Cushion

ICEROSS DERMO® CUSHION LINER





Icross Dermo Cushion liner is made with the softest silicone that actively cares for the skin.

This liner provides superb softness, gentle skin contact and Össur's unique Active Skin Care ingredients. The combination of DermoGel® silicone with an ultra-strong and elastic Supplex® fabric cover offers excellent durability and an intimate fit.

Össur recommends that Icross Dermo Cushion Liner is used in conjunction with Icross Sleeve and Icelock® 500 Series.

USER INFORMATION	
Amputation Level:	Transibial
Impact Level:	Low to High
LINER INFORMATION	
Size Standard:	16, 18, 20, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 40*
Profile:	3mm or 6mm

*Size 40 is only available in 3mm.



SUSPENSION METHOD



FEATURES



PART NUMBER INFORMATION - ICEROSS DERM O CUSHION LINER

PART#	PROFILE	SUSPENSION METHOD
I-CL43XX	3mm	Cushion
I-CL46XX	6mm	Cushion