


ICEROSS DERMO® JUNIOR LOCKING LINER WITH WAVE



Iceross Dermo Junior Locking Liner with Wave offers superb comfort and impact protection for children with sensitive skin engaging in medium impact daily activities. Active Skin Care ingredients are built into the silicone material.

Children can decorate their liners and personalize them. There are 4 different colored permanent markers included with the liner.

Össur recommends that Iceross Dermo Junior Locking Liner with Wave is used in conjunction with the Icelock® 700 Series.



USER INFORMATION

Amputation Level:	Transtibial application for children
Impact Level:	Low to High

LINER INFORMATION

Size Standard:	12, 14, 16, 18.
Profile:	3mm
Matrix:	10mm

SUSPENSION METHOD



LOCKING

FEATURES



MATRIX



FABRIC COVER



ACTIVE SKIN CARE



WAVE



SILKEN INNER SURFACE

PART NUMBER INFORMATION - ICEROSS DERMO JUNIOR WAVE LINER

PART#	PROFILE	SUSPENSION METHOD	MATRIX
I-8302XX	3mm	Locking	10cm

Note: Sizes 12 and 14 use M6 thread. Sizes 16 and 18 offer interchangeable distal pin attachment for use with M6 or M8

