

ICELOCK® 700 SERIES




Designed to be used in conjunction with the Iceross Upper X and Iceross Junior liner systems.

- Extremely low build height and compact design
- Light-weight
- Compatible for single or double lamination fabrication method
- 45kg (99lbs) weight limits on components
- Locking mechanisms are replaceable from inside of the socket
- All impact levels



DIAGRAM



Double lamination
Icelock Ratchet 721



Single lamination
Icelock Stainless Steel Pyramid 773



Single lamination
Icelock Alloy 4-hole 774

PART#	DESCRIPTION	RATINGS		WEIGHT	
L-721000	Icelock 721 Ratchet, Small (No Adapter)*	13mm	1/2"	20g*	0.71oz

*Use with L-780000 Fabrication Tooling Kit

PART#	DESCRIPTION
L-780000	Icelock fabric tooling kit

SPECIALLY DESIGNED PYRAMIDS FOR ICELOCK 700

PART#	DESCRIPTION	RATINGS		BUILD HEIGHT	
L-773000	Icelock 773 Stainless Pyramid	45kg	99lbs	13mm	1/2"

PART#	DESCRIPTION	RATINGS		BUILD HEIGHT	
L-774000	Icelock 774 Aluminum 4-Hole	45kg	99lbs	13mm	1/2"

RELEASE KEY

PART#	DESCRIPTION
L-792016	A tool to tighten and release the release unit (sold separately)



ICELOCK® 700 SERIES

COMPONENTS FOR RATCHET LOCK 721

PART#	DESCRIPTION
L-792150	Fabrication Ring, small
L-792020	Ratchet Lock without Adapter, small
L-792110	Funnel, small
L-792027	Release Unit, small
L-792031	Pin Adapter
L-792100	Attachment Pin Ratchet, 32mm
L-792105	Attachment Pin Ratchet, 56mm
L-793112	Silicone Seal, small

COMPONENTS FOR FABRICATION (HARDENED STAINLESS STEEL)

PART#	DESCRIPTION
L-792119	Steel Guide, small
L-792010	Lock Dummy with o-ring, small
L-792120	Release Dummy, small
L-792130	Pin Dummy, small
L-793110	Pin Dummy Bolt, small
L-692015	Hex Driver with o-ring
L-792017	Cross Driver with o-ring, small