


ICEROSS® ACTIVA WITH TIBIAGUARD™ CUSHION LINER



The supreme cushion liner for active users with optimal cushioning, comfort and maximum stability.

The Iceross Activa with TibiaGuard is manufactured with revolutionary techniques that combine different types of silicone to meet the challenging requirements of highly active transtibial users.

The integrated TibiaGuard provides optimal cushioning while a thinner, stiffer posterior portion is designed to enhance maximum stability.



USER INFORMATION	
Amputation Level:	Transtibial
Impact Level:	Moderate to Extreme
LINER INFORMATION	
Size Standard:	18,20,22,23.5,25,26.5,28,30,32,34.
Profile:	2mm posterior, 6mm anterior

SUSPENSION METHOD

FEATURES



CUSHION



FABRIC COVER



ACTIVE SKIN CARE



WAVE



DOUBLE DURO



SILKEN INNER SURFACE



TIBIAGUARD

PART NUMBER INFORMATION - ICEROSS ACTIVA WITH TIBIAGUARD

PART#	PROFILE	SUSPENSION METHOD	ITEMS
I-CT62XX	2mm / 6mm	Cushion	N/A