

CTI®

The proven technology of the CTi knee brace is available with off-the-shelf convenience, providing adjustability and an excellent out-of-the-box fit while maintaining the level of support and protection that you expect from a CTi knee brace.

Outfitted with anatomically-correct Accutrac® hinges and flexible subshells and d-rings, the rock-solid CTi carbon frame provides a rigid exoskeleton that stabilizes the knee joint. Breathable liners coated with bio-inert Össur Sensil® Silicone ensure CTi OTS remains properly positioned on the leg, for optimal functional benefit.



INDICATIONS FOR USE

- For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities

FEATURES AND BENEFITS

- Anatomically-correct Accutrac hinges with extension stops to glide with the knee
- Flexible subshell and Adjustable Tibial Strap capture the tibia and allow adjustments in hinge depth for a more customized fit
- Breathable liners coated with bio-inert Össur Sensil Silicone to reduce migration and skin irritation
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Non-corrosive materials, ideal for water sports
- Available in 2 models: Standard model for moderate to high impact levels, Pro Sport model for the highest impact levels

REQUIRED MEASUREMENTS

- Caliper measurement: Medio-laterally at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella



CTI®

ORDERING INFORMATION

CTI

Size	Caliper	Circumference	Standard Left	Standard Right	Pro Sport Left	Pro Sport Right
Small	8.9 – 10.2 cm	28 – 34 cm	B-238600112	B-238500112	B-239600112	B-239500112
Medium	10.2 - 11.4 cm	34 – 38 cm	B-238600113	B-238500113	B-239600113	B-239500113
Large	11.4 – 12.7 cm	38 – 41 cm	B-238600114	B-238500114	B-239600114	B-239500114
X-Large	12.7 – 14.0 cm	41 – 44 cm	B-238600115	B-238500115	B-239600115	B-239500115
XX-Large	14.0 – 15.2 cm	44 – 48 cm	B-238600116	B-238500116	B-239600116	B-239500116

Note: M/L measurements are the preferred and most accurate form of measurement for the CTI. Circumferential measurements are approximate.

ACCESSORIES

LINERS

Size	Silicone tibia cuff	Silicone calf	Silicone thigh	Lycra thigh
Medium	B-304030030	B-425580002	B-425571002	B-425570002
Large		B-425580004	B-425571004	B-425570004
X-Large, XX-Large		B-425580005	B-425571004	B-425570004

STRAPS KITS

Size	Part#
Small, Medium	B-701140103
Large, X-Large, XX-Large	B-701140104

PCL KITS

Size	Part#
Small, Medium	B-705000132
Large, X-Large, XX-Large	B-705000235

STRAP PADS

Size	Part#
Small	B-304010057
Medium, Large	B-304010058
X-Large, XX-Large	B-304010059

AMS (ANTI-MIGRATION SYSTEM) WRAPS

Description	Part#
Silicone AMS Wrap Small, Medium	B-715260002
Silicone AMS Wrap Large, X-Large, XX-Large	B-715260004
Neoprene AMS Wrap Small, Medium	B-410000001
Neoprene AMS Wrap Large, X-Large	B-410000002

CTI®

ACCESSORIES

MX GEAR

Description	Part#
Patella Cup	B-711130011
Gear Guards	B-705000003
MX Kit S (patella cup, gear guards, undersleeve, stickers)	B-711140002
MX Kit M (patella cup, gear guards, undersleeve, stickers)	B-711140003
MX Kit L (patella cup, gear guards, undersleeve, stickers)	B-711140004
MX Kit XL (patella cup, gear guards, undersleeve, stickers)	B-711140005

SLEEVES

Circumference 6" (15 cm) above mid patella	Softsleeve	Undersleeve	Oversleeve	Sportsleeve
<15" (38 cm)	20643	20251	20051	20682
15 – 18" (38 – 45 cm)	20644	20252	20052	20683
18 – 19.5" (45 – 49.5 cm)	20644	20253	20053	20684
19.5 – 21.5" (49.5 – 54.5 cm)	20645	20254	20054	20685
21.5 – 24.5" (54.5 – 62 cm)	20645	20255	20055	N/A
24.5 – 28" (62 – 71 cm)	20646	20256	20056	N/A
28 – 34" (71 – 86 cm)	20647	20257	20057	N/A

OTHER ACCESSORIES

Description	Part#
CTi D-Ring Kit Standard 2x 1" and 2x 1.5" d-rings	B-704132600
CTi D-Ring Kit 2" 2x 2" d-rings	B-704132602
CTi Condyle Pad Kit	B-704130150
Extension Stop Kit	22130
Flexion Stop Kit	19254