

SHIPPING INFORMATION

SHIP VIA: GROUND NEXT DAY 2ND DAY

CALL TO CONFIRM ORDER?: YES NO

DATE: _____ ACCT.#: _____ PURCHASE ORDER #: _____

REQUIRED PATIENT INFORMATION

LAST NAME: _____ FIRST NAME: _____

HEIGHT: _____ ft. _____ in. WEIGHT: _____ lbs. FOOT SIZE: _____ cm.

ACTIVITY LEVEL: K1 K2 K3 K4 CATEGORY (1-9): _____ CLEARANCE: _____ FOOT SHELL COLOR: BEIGE BROWN

AMPUTATION TYPE: LEFT RIGHT BILATERAL LEVEL OF AMPUTATION: HD KD SYMES BK AK

IMPACT LEVEL: HIGH: Daily activities involving rigorous and repetitive actions at high levels of impact. MEDIUM: Daily activities involving normal and repetitive actions at moderate levels of impact. LOW: Daily activities involving limited and steady actions at low levels of impact.

BILLING ADDRESS

Contact: _____
Company: _____
Address: _____
City: _____ State/Zip: _____
Phone: _____ Fax: _____

SHIPPING ADDRESS (if different)

Contact: _____
Company: _____
Address: _____
City: _____ State/Zip: _____
Phone: _____ Fax: _____

	PRODUCT NAME	AVAILABLE ADAPTER OPTIONS					QTY	NOTES (Ex: Heel Height, Solid Toe, etc.)
		MALE PYRAMID	FEMALE PYRAMID	TUBE (Short or Long)	4-HOLE	T ADAPTER		
IMPACT SOLUTIONS	Re-Flex Rotate - Standard							
	LP Rotate							
	Re-Flex Shock - Standard							
DYNAMIC SOLUTIONS	Pro-Flex® LP							
	Pro-Flex® XC							
	PROPRIO FOOT®							
	LP Vari-Flex							
	Vari-Flex®							
	Unity Pump Module - Size 22-24							
	Unity Pump Module - Size 25-27							
BALANCE SOLUTIONS	Unity Pump Module - Size 28-30							
	K2 Sensation®							
	Balance™ Foot J							
	Flex-Foot® Assure							
	Unity Pump Module - Size 23-26							
ICELOCK 500	Unity Pump Module - Size 27-28							
	Valve							<input type="checkbox"/> TF <input type="checkbox"/> TT
	544 Unity Plate							

	PRODUCT NAME	PROFILE	SIZES	QTY	NOTES
LINERS	Iceross® Seal-In® V	3MM	18,20,22,23.5,25,26.5,28,30,32,34,36		
	Iceross® Seal-In® V HP	3MM	18,20,22,23.5,25,26.5,28,30,32,34,36		
	Iceross® Seal-In® X	<input type="checkbox"/> 3MM <input type="checkbox"/> 6MM	Standard: 16,18,20,22,23.5,25,26.5,28,30,32,34, 36,40 Conical: 23.5,25,26.5,28,30,32,34,36,38,40 Seal Ring: 20,22,24,26,28,30,32,35,38,41		
	Iceross® Seal-In® X TF	3MM	25,26.5,28,30,32,34,36,38,40,45,50,55 Seal Ring: 26,28,30,32,35,38,41,44,47,51,55,60,65,70		

* The 3/4" heel and carbon tube options are not yet available for Flex-Foot feet when combined with Unity Sleeveless Vacuum.

Measure for correct liner size 4 cm up from the distal end of the limb. Do not downsize for this seal-in liner version. Order the size closest to your measurement.

If you require a custom liner, please contact Evolution Industries, an Össur company, for two other validated Seal-In options for our Unity System.

888-839-6213 (US & Canada) or cs@evoi.com.