



REBOUND® PCL

The Rebound PCL is the first brace designed to apply a physiologically correct, dynamic force, optimum for rehabilitation of PCL ligament. Insufficient healing after a PCL rupture can lead to significant pathologies, such as, chronic knee pain or knee OA.¹

The Dynamic Tension System (DTS)¹ is designed to increase the load on the tibia as the knee goes through flexion providing an Anterior Drawer force necessary to decrease load on the PCL to allow it to heal. The Rebound PCL allows for specific load adaptation² (4, 7, or 10 lbs) according to the patient's individual anatomy and rehabilitation related requirements. The specially designed Anit-migration Straps (AMS)³ reduces migration and provides comfort to the patient while ensuring that the brace stays in the correct position during wear.



INDICATIONS

- Designed to help optimize the healing process following posterior cruciate ligament (PCL) ruptures, both during functional (non-surgical) and post-surgical rehabilitation.

PART#	DESCRIPTION
B-124501154	Rebound PCL Brace Custom Right
B-124501155	Rebound PCL Brace Custom Left
PCL00010	Extension Stop Kit
PCL00020	Flexion Stop Kit
PCL00030	Shear Knob Kit
PCL00061	S-L Strap Kit
PCL00062	XL-XXL Strap Kit



Nothing measures like SmartMeasure.™

An app for iPhone® & iPad® designed to measure custom knee braces, SmartMeasure is Össur's most intuitive, efficient and accurate way to measure for any Össur custom knee brace.

Apple, the Apple logo, iPhone and iPad are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.



FUNCTIONAL HEALING™ – WHERE MOBILITY MEANS RECOVERY

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimize healing while maintaining function and mobility.