



UNLOADER® HIP

The Unloader Hip is designed to optimize Hip OA treatment in order to improve mobility and decrease pain. It's unique Dynamic Rotation Straps™ is clinically proven to decrease the internal rotation and adduction moment. It also delivers proprioceptive support during the gait cycle by providing compression with a proven pulley system design.




FEATURES AND BENEFITS

- 1 The Dynamic Rotation Straps allows individual adjustment of the external rotation force, to achieve optimum pain relief for each patient
- 2 With its simple-to-tighten design, the proven pulley system provides compression and proprioceptive support
- 3 The trochanter pad helps relieve pressure on the greater trochanter region

INDICATIONS

- For mild to moderate hip osteoarthritis (uni-or-bilateral)



SIZE	MEASUREMENTS ¹	WOMEN'S	MEN'S
Small	31.5-35" (80-88cm)	B-465350030	B-465350031
Medium	35-38" (88-96cm)	B-465350032	B-465350033
Large	38-41" (96-104cm)	B-465350034	B-465350035
XLarge	41-45" (104-114cm)	B-465350036	B-465350037
XXLarge	45-49" (114-124cm)	B-465350038	B-465350039

ACCESSORIES DESCRIPTION

DESCRIPTION	PART#
B-465350064	Rotation Strap Kit
B-465350065	Washbag Kit
B-465350066	Trochanter Pad Kit

SIZE	THIGH STRAP KIT	BILATERAL STRAP KIT
Small	B-465350067	B-465350072
Medium	B-465350068	B-465350073
Large	B-465350069	B-465350074
XLarge	B-465350070	B-465350075
XXLarge	B-465350071	B-465350076

SIZE	WOMEN'S TROUSER	MEN'S TROUSER
Small	B-465350077	B-465350078
Medium	B-465350079	B-465350080
Large	B-465350081	B-465350082
XLarge	B-465350083	B-465350084
XXLarge	B-465350085	B-465350086

¹ Measure the hip circumference at the anterior superior iliac spine (hip circumference). See drawing.