

## PRO-FLEX®



Combining an incredible 27° ankle motion and significantly greater ankle power than a conventional carbon fiber foot, Pro-Flex has been developed specifically to enhance protection of the sound side. It closely mimics regular biomechanics, delivering a proven 11% reduction in load on the contralateral limb.

- 27° ankle motion
- Helps reduce load on contralateral limb by 11%
- Exceptional energy return
- Innovative pivot mechanism
- Aesthetically pleasing, anatomical footcover

Össur recommends using Pro-Flex in conjunction with the sleeveless vacuum suspension system Unity®, which adds minimal weight and no extra build height, and does not interfere with Pro-Flex dynamics. Unity is designed to provide secure suspension in combination with Iceross Seal-In® technology, while stabilizing volume fluctuations and maintaining comfortable socket fit throughout the day.



User Information	
Amputation Level:	Transfemoral and Transtibial
Impact Level:	Low to Moderate
Maximum Patient Weight:	125 kg (275 lbs)

Unity System Information	
Maximum Vacuum:	-23inHg (-600mmHg)
Normal Operating Pressure:	-16-20 inHg
Added System Weight:	75g (0,17 lbs) w/ Pump, Socket Valve and Tube


Foot information	
Categories:	1-7
Sizes:	22-30
Weight of Foot: (Size 27, Cat. 5)	920 g (2,02 lbs) w/ Foot Cover
Build Height: (Size 27)	155 mm (6 1/8")
Heel Height:	10mm (3/8")
Adapter Options:	Male Pyramid



# PRO-FLEX®

## ORDERING INFORMATION


### FOR COMPLETE FOOT



Please select the appropriate Unity Pump Module kit according to foot module size when placing an order and Unity will come preassembled on the foot. For extra convenience and system compatibility, Unity Valve Kits are highly recommended as a means of connecting to the socket.

For liner selection charts and item numbers please refer to the appropriate Icross Seal-In page.

### FOR COMPLETE FOOT WITH UNITY




## CATEGORY SELECTION CHART

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-125
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-275
<b>Low Impact Category</b>								
Sizes 22-24	1	1	2	3	4	5	N/A	N/A
Sizes 25-27	2	2	2	3	4	5	6	7
Sizes 28-30	3	3	3	3	4	5	6	7
<b>Moderate Impact Category</b>								
Sizes 22-24	1	2	3	4	5	N/A	N/A	N/A
Sizes 25-27	2	2	3	4	5	6	7	N/A
Sizes 28-30	3	3	3	4	5	6	7	N/A

# PRO-FLEX®



## CLEARANCE



\*Note: Pro-Flex with Unity has the same build height as Pro-Flex.

## KITS

### PRO-FLEX



Includes:

- Foot Module
- Removable Foot Cover\* with a black Spectra Sock

\*For Pro-Flex with brown foot cover add a "BR" suffix to the part#.

# PRO-FLEX®

## UNITY

### UNITY KITS

Picture	Part#	Description	Applicable Foot Kits
1	UPF0001	PFP Unity Pump S22-30	Pro-Flex Size 22-30
2	UTT0001	Unity TT Valve Kit	All
3	UTF0001	Unity TF Valve Kit	All
4	UTF0006	Unity TF Tube Kit	All
5	L-5446XX	Icelock 544 Kit	All




### UNITY SPARE COMPONENTS

Part#	Description
UPF00011	Unity PFP Components 22-30

Includes:


- Bumper
- Membrane
- Mounting screws
- Loctite



### FOOT KITS APPLICABLE WITH UNITY


#### PRO-FLEX

PART#



Includes:

- Foot Module
- Removable Foot Cover\* with a black Spectra Sock




\*For Pro-Flex with brown foot cover add a "BR" suffix to the part# Note: Unity will come preassembled on the foot if ordered with PFPU

# PRO-FLEX®

## COSMETIC FINISHING

### FOOT COVER KIT, BEIGE AND BROWN\*

PART#



Includes:

- Footcover
- Footcover Attachment
- Black Spectra Sock

\*Note: Please specify when ordering brown covers, add "BR" to the part#.