

RHEO KNEE® XC

UNITY
COMPATIBLE

INTRODUCTION

From rehabilitation, to the home and office, to the gym and the trail, RHEO KNEE XC confidently goes anywhere. Whether rehabilitating on a spin cycle, navigating tight spaces, strolling around town, or hiking a trail in the rain, RHEO KNEE XC provides the perfect balance of stability and dynamics.

KEY FEATURES

- Weatherproof design allows exposure to splashing of water.
- Magnetorheologic technology enables an instant response so that users never have to wait for the knee to catch up with them.
- Smart gait detection, including kinematic sensor technology, ensures stability and dynamic response in every situation.
- Advanced actuator and resistance control ensure best possible resistance, e.g. more support in stair descent and minimum effort needed in level ground gait.
- Advanced extension assist mechanism.
- Effortless swing initiation enables a smoother gait, even in crowds and confined spaces.
- Allows high impact activities such as running and jogging.
- Extension lock mechanism, to lock the knee in full extension and increase safety and comfort under specific circumstances.
- Efficient cycling and secure stops with Automatic Cycling Recognition.
- Allows step-over-step stair ascent.
- ÖSSUR LOGIC app offers two levels of experience one for the user and another for the practitioner.

RECOMMENDED WITH



Seal-In® X TF liner



Pro-Flex®



Re-Flex® Shock



Re-Flex® Rotate



Pro-Flex® XC

*Note. Pro-Flex can not be used for high impact activities.



ÖSSUR KNEES
AND LEGS



RHEO KNEE® XC




User Information		
Amputation Level:	Transfemoral and Knee Disarticulation	
Impact Level:	Moderate	High
Maximum Patient Weight:	136kg (300lbs)	110 kg (243 lbs)


Knee Information	
Net weight:	Approx. 1.61 kg (3.56 lbs)
Build Height:	236mm (9 1/4")
Knee Flexion:	120°
Aluminum Frame Construction	
Weatherproof	

CLEARANCE

RHEO KNEE XC



RHEO KNEE XC WITH THREADED TOP ADAPTER (A - 145310)



RHEO KNEE® XC

KITS

RHEO KNEE XC KIT

Part#	Description
RKNXC0002	2 Years Warranty (No Service Check)
RKNXC0003	3 Years Warranty (No Service Check)
RKNXC0004	4 Years Warranty (Incl. Service Check at 40 Months)
RKNXC0005	5 Years Warranty (Incl. Service Check at 40 Months)

Includes:

- RHEO KNEE XC
- RHEO KNEE Power Supply w/Plug
- Warranty Certificate
- User Manual



RHEO KNEE XC KIT WITH THREADED TOP ADAPTER:

RHEO KNEE XC with threaded top adapter is the solution for knee disarticulation amputees and those with long residual limbs.

To order RHEO KNEE XC with threaded top adapter, order the RHEO KNEE XC Kit and add Össur's Female Pyramid Insert to the order.

Part#	Description	Height		Weight	
A-145310	Female Pyramid Insert / Threaded Top Adapter	9mm	3/8"	33g	1.16oz

RHEO KNEE ACCESSORIES

Part#	Description
RKNLPCWB1	RHEOLOGIC Workbench Software for PC use
N/A	ÖSSUR LOGIC for iOS devices is available on the App Store

Part#	Description
BIC00110	Bionic Protector for RHEO KNEE



RHEO KNEE® XC

COMPONENTS

POWER SUPPLY

Part#	Description
PRX00421	Power Supply w/Plug



EXTENDED WARRANTY

RHEO KNEE XC KIT EXTENDED WARRANTY*

Part#	Description
RKNXCGE1YR	Additional 1 Year Warranty
RKNXCGE2YR	Additional 2 Year Warranty (Incl. Service Check at 40 Months)
RKNXCGE3YR	Additional 3 Year Warranty (Incl. Service Check at 40 Months)
RKNXCSC20M	Voluntary 20 Month Service Check
RKNSCAW	Voluntary After Warranty Service Check

* Maximum total warranty period is 5 years, must be purchased within 20 months after original purchase date.

UNITY

UNITY KITS

Part#	Description	Applicable Foot Kits
UTF0001	Unity TF Valve Kit	All

