



CHEETAH® KNEE

Lean and mean. Like you.

As an athlete, you need to be as efficient as possible. Energy can't be wasted and every ounce matters. That's why, in creating the Cheetah Knee, we stripped our rugged Total Knee 2100 to a lean 22 ounce core, removing everything but the essentials. Then, we redesigned the frame and linkages, and optimized the hydraulics for faster flexion/extension. So, whatever your finish line, you get there as quickly as possible.



"The lightest, fastest knee I've ever run on. Cheetah Knee is a game changer."

Sarah Reinertsen



CHEETAH® KNEE

- For jogging and distance running
- Short knee center to top accommodates longer residual limbs
- 4-bar geometry with stable stance control for controlled deceleration
- The perfect match for Flex-Run™ (see below)
- Enhanced hydraulics for faster swing speed
- Adjustable 3-phase hydraulic swing control
- Extension bumper for softer terminal impact

RATINGS

- Rated for high impact
- Weight limit of 220lbs (100kg)
- Light weight, only 617g
- 25mm knee center to proximal threaded top, ideal for long residual limb TF users
- Build height 108mm
- Flexion: 110°

OPTIONS

- CKN10000 Cheetah Knee

18-Month Warranty

SUGGESTED BILLING CODE*
L5613, 5930



* The responsibility for accurate coding lies solely with the provider treating the patient. Össur assumes no responsibility or liability for the provider's coding decisions. Össur's coding suggestions rest on its best judgment and are subject to revision.



Designed to empower athletes to fulfill their potential in sporting events and activities.

FLEX-RUN WITH NIKE® SOLE

For jogging and competitive distance running on or off road. Includes removable Nike Sole.



FOLLOW ÖSSUR ON



USA (800)233-6263
CANADA (800)663-5982
WWW.OSSUR.COM

