

ÖSSUR®

FORM FIT®

THE CHOICE OF EXPERTS



Form Fit Knee Brace

Sometimes, our bodies need assistance. For decades, specialists have trusted Form Fit supports, a complete range of quality products by Össur. A global leader in innovative orthopaedic products, Össur's biomechanical expertise has helped millions live Life Without Limitations.

That's why Form Fit is the Choice of Experts.



FORMFIT®

THE CHOICE OF EXPERTS



FORM FIT® KNEE BRACE

The Form Fit Knee brace is for patients suffering from mild to moderate strains, sprains or ligament tears (ACL/PCL/MCL/LCL). It provides excellent stability, protection and compression during injury recovery. Made from strong, breathable materials, it offers comfort while remaining highly durable. The open popliteal area reduces material bulging and potential irritation, while the malleable rigid hinge arms provide an intimate fit. Available in a wrap or sleeve version. ROM available in 12" or 16" length.

Features and Benefits

- Breathable material
- Easy pull-on
- Malleable Hinges
- Universal buttress
- Open popliteal

Indications

- Mild to moderate ACL and PCL tears
- MCL and LCL tears
- Combined mild knee sprains and strains
- Instabilities of the knee

Size	Circumference	Rom		Non-Rom	
		Sleeve	Wrap	Sleeve	Wrap
XSmall	13-15.5" (33-39.4cm)	508X52	507X52	504252	503252
Small	15.5-18.5" (39.4-47cm)	508X53	507X53	504253	503253
Medium	18.5-21" (47-53.3cm)	508X55	507X55	504255	503255
Large	21-23.5" (53.3-59.7cm)	508X57	507X57	504257	503257
XLarge	23.5-26.5" (59.7-67.3cm)	508X58	507X58	504258	503258
2XLarge	26.5-29.5" (67.3-74.9cm)	508X59	507X59	504259	503259
3XLarge	29.5-32" (75-81cm)	—	507X60	—	503260

X = 3 (Long) or 2 (Short)

*Thigh Circumference: 6" (15cm) above mid-patella

ROM Hinge Flexion & Extension Stops

Extension Stops: 10°, 20°, 30°, 40°, lockout (full extension).

Flexion Stops: 45°, 60°, 75°, 90°.

FOLLOW ÖSSUR ON



TEL (800) 233-6263
CANADA (800) 663-5982
WWW.OSSUR.COM

