

UNLOADER ONE® LITE

Unloader One Lite is a slim and low profile knee brace that relieves pain from degenerative meniscal tears and knee OA by means of the clinically proven Dynamic Force System that provides uni-compartmental unloading of the knee joint.

With an integrated Lycra sleeve and shells with overmold, the brace is an optimal combination of comfort and sleek look.



INDICATIONS FOR USE

- Mild to moderate unicompartmental osteoarthritis
- Degenerative meniscus
- Other knee conditions that may benefit from unicompartmental unloading

FEATURES AND BENEFITS

- 3-Point Leverage System provides clinically relevant pain relief of affected compartment
- Single SmartDosing dial allows fine-tuned pain management
- Helps keep patients mobile during daily activities
- Discreet hinge and integrated sleeve ensure a sleek look on leg and snug fit underneath clothing
- Lightweight and easy to use

REQUIRED MEASUREMENTS

- Circumference measurement: 6"/15 cm above and below mid-patella



UNLOADER ONE® LITE

ORDERING INFORMATION

UNLOADER ONE LITE

| Size | Circumference 6" (15 cm) above mid patella | Circumference 6" (15 cm) below mid-patella | Left medial | Right medial | Left lateral | Right lateral |
|---------|--|--|-------------|--------------|--------------|---------------|
| X-Small | 12 - 15" (30 - 38 cm) | 10 - 12" (25 - 30 cm) | UOL23000 | UOL13000 | UOL24000 | UOL14000 |
| Small | 15 - 17.75" (38 - 45 cm) | 12 - 13.75" (30 - 35 cm) | UOL23010 | UOL13010 | UOL24010 | UOL14010 |
| Medium | 17.75 - 22" (45 - 56 cm) | 13.75 - 16.5" (35 - 42 cm) | UOL23020 | UOL13020 | UOL24020 | UOL14020 |
| Large | 22 - 26" (56 - 66 cm) | 16.5 - 20" (42 - 51 cm) | UOL23030 | UOL13030 | UOL24030 | UOL14030 |
| X-Large | 26 - 29" (66 - 73 cm) | 20 - 24" (51 - 61 cm) | UOL23040 | UOL13040 | UOL24040 | UOL14040 |

Hinge is placed on the side of the affected compartment