

REBOUND® POST OP KNEE

The simple, versatile post-operative knee bracing solution that is easy for everyone. Whether adjusting ROM, setting the drop lock, or altering the brace length, you'll appreciate the simplicity of the Rebound Post-op knee brace.



INDICATIONS FOR USE

- ACL, PCL, MCL, LCL repairs/injuries
- Meniscal repairs
- Tibial plateau fractures
- Patellar tendon repairs
- Osteochondral repairs
- Condylar fractures
- Sprains/strains of the knee
- HTOs (High Tibial Osteotomy)

FEATURES AND BENEFITS

- Tool-free design and multi-purpose Smart-Fit™ packaging
- Available in universal, left or right contour and pediatric
- Telescoping frame from 18.5"-27" (Adult) and 15"-21" (Pediatric)
- OneTouch™ easy to adjust ROM hinge with lock-out capability
- Color-coding and numbered frame for patient education
- Quick-fit, gripping buckles for easy application

REQUIRED MEASUREMENTS

- Circumference: 6" (15 cm) above mid-patella

ORDERING INFORMATION

REBOUND POST-OP KNEE

Description	Thigh circumference	Cool	Full foam
Pediatric	12 – 25" (30,5 cm – 63,5 cm)	B-8381470	
Universal, Left contoured	16 – 32" (41 cm – 81 cm)	B-8381680	B-8382680
Universal, Right contoured	16 – 32" (41 cm – 81 cm)	B-8381580	B-8382580
Universal	16 – 32" (41 cm – 81 cm)	B-8381480	B-8382480
Universal XL	16 – 38" (41 cm – 96,5 cm)	B-8381490	B-8382490

ACCESSORIES

Description	Part#
ROM stop locks	B-8380005
Alligator velcro kit	B-8380006
Adhesive velcro kit pediatric	B-8380008
Adhesive velcro kit universal / universal XL	B-8380007

REPLACEMENT SOFTGOODS

Description	Cool	Full foam
Pediatric	B-8381072	
Universal*	B-8381082	B-8382082
Universal XL	B-8381092	B-8382092

*For universal, left contoured, right contoured