

ÖSSUR FORMFIT® PRO KNEE

Össur Formfit® Pro Knee is a compressive knee sleeve that features a SmartPad™ buttress for targeted support and cushioning. Our patented design realigns the patella, providing pain relief from conditions such as Patella Tendinitis (Jumper’s Knee) and Chondromalacia Patella (Runners’ Knee). The dynamic, 3-D MotionTech™ knit is supportive, meets certified medical grade compression requirements and has been proven to be superior to the leading brands at moisture wicking and breathability, keeping your knee cool during wear. Light weight, removeable spiral stays ensure optimum shape and fit, while distinct patella graphics make donning intuitive.



PROKNEE



MOTIONTECH



COOLVENT



STAYBILIZER



SMARTPAD



FOLDLOCK

INDICATIONS FOR USE

- For knee conditions that may benefit from compression, such as;
- Patella tendinitis
- Chondromalacia Patella
- Mild sprains and strains
- Knee pain
- Feeling of instability

CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your knee cool during wear
- Certified medical grade compression (20-30 mmHg)
- Lightweight CoolVent™ knit in popliteal area for added breathability and comfort
- Patented SmartPad™ for targeted support and cushioning
- Patella graphics for intuitive placement
- Patented FoldLock™ system for easier and more secure application
- Lightweight, removable spiral steel stays

REQUIRED MEASUREMENTS

- Calf circumference 15 cm (6") below mid-patella
- Thigh circumference 15 cm (6") above mid-patella

Size	Calf		Thigh	
XS	26 – 31 cm	10 ¼ – 12 ¼ inch	36 – 41 cm	14 ¼ – 16 ¼ inch
S	31 – 36 cm	12 ¼ – 14 ¼ inch	41 – 46 cm	16 ¼ – 18 inch
M	36 – 41 cm	14 ¼ – 16 ¼ inch	46 – 51 cm	18 – 20 inch
L	41 – 46 cm	16 ¼ – 18 inch	51 – 56 cm	20 – 22 inch
XL	46 – 50 cm	18 – 19 ¾ inch	56 – 60 cm	22 – 23 ½ inch

KNEE

ÖSSUR FORMFIT® PRO KNEE

ORDERING INFORMATION

ÖSSUR FORMFIT® PRO KNEE BLACK

Size	Part#
XS	FP100001
S	FP100002
M	FP100003
L	FP100004
XL	FP100005



ÖSSUR FORMFIT® PRO KNEE BLUE

Size	Part#
XS	FP101001
S	FP101002
M	FP101003
L	FP101004
XL	FP101005

