

## DH OFFLOADING WALKER®

### INDICATIONS FOR USE

- Treatment of plantar ulcers

### FEATURES AND BENEFITS

- Shear force reduction by stabilizing foot and ankle at 90°
- Wide low profile rocker bottom provides a smooth roll-through and increased stability for altered gaits
- The softgood's patented inner bootie prevents foot from migrating on the insole
- Designed around Össur's unique Active Offloading Insole
- Provides a safe environment for the simultaneous use of wound care material
- Easy to apply and modify

### REQUIRED MEASUREMENTS

- Shoe size



### ORDERING INFORMATION

#### DH OFFLOADING WALKER

Size	Shoe Size US (Men)	Shoe Size US (Women)	Shoe Size EU	Item No
X-Small	2-4	3.5-5.5	33-35.5	DH0200BLK
Small	4.5-7	6-8	36-39	DH0400BLK
Medium	7.5-10.5	8.5-11.5	39.5-44	DH0600BLK
Large	10.5-12.5	11.5-13.5	44.5-46.5	DH0800BLK
X-Large	12.5+	13.5+	47-50	DH1000BLK

### ACCESSORIES

#### REPLACEMENT PARTS

Size	Softgoods	Insoles
X-Small		11742C
Small	10958	PP04C
Medium	10959	PP06C
Large	10960	11740C
X-Large		11741C

## DH OFFLOADING POST-OP STYLE SHOE

### INDICATIONS FOR USE

- Plantar wounds

### FEATURES AND BENEFITS

- Designed to accommodate either the left or the right foot
- Designed around Össur's patented Active Offloading Insole
- Designed to accommodate bulky dressings
- Super-soft nylon upper with reinforced heel counter provides extra comfort and stability
- Patented removable toe piece aides in wound accessibility for daily wound site inspection

### REQUIRED MEASUREMENTS

- Shoe size



### ORDERING INFORMATION

#### DH POST-OP SHOE

Size	Shoe Size US (Men)	Shoe Size US (Women)	Shoe Size EU	Item No
Small	4-6	6.5-8.5	36-39	10340
Medium	6.5-8.5	9-11	39.5-44	10341
Large	9-11	12+	44.5-46.5	10342
X-Large	11.5-14	13+	47-50	10343

### ACCESSORIES

#### REPLACEMENT PARTS

Size	Insoles
Small	12365
Medium	12366
Large	12367
X-Large	12368