



ÖSSUR[®]
LIFE WITHOUT LIMITATIONS

STAY ALIGNED WITH THE PRO-FLEX LP ALIGN

Regardless of whether you want smart shoes, comfortable sneakers, or simply to kick off your shoes after a long day and walk barefoot at home, it's important to stay aligned.

DRESS FOR THE OCCASION, WITHOUT COMPROMISE

The new Pro-Flex LP Align tackles the challenge of changing shoes head on by empowering you to align your prosthesis yourself, ensuring that you don't have to compromise your posture or gait because of different footwear. Pro-Flex LP Align allows you to align your prosthesis with the shoes you wear with the simple touch of a button, up to 7 cm heel height!

MORE STABILITY AND LESS FATIGUE TO BOOST YOUR CONFIDENCE

Pro-Flex LP Align not only offers easy to use adjustable alignment but also has good energy return and helps to reduce fatigue, whilst boosting stability and confidence. The anatomically-shaped foot blade with a sandal toe and full effective toe length design offer enhanced stability and better roll over, and contribute to a more symmetrical gait. The slight curvature at the toe also ensures a good fit at different heel-height settings.

ANATOMICAL SHAPED COVER FOR OPTIMAL FUNCTION

Lightweight and washable, the foot cover comes with a gripping sole for barefoot stability on wet or slick surfaces, plus a 'sandal toe' for use with flip-flops or sandals. A calf cosmesis made from the same light weight and durable material as the foot cover is offered for sizes 22-26 (cat 1-5), specially to tolerate the alignment.

HOW TO STAY ALIGNED


To prevent mal-alignment and its resulting problems, Össur has created a video tutorial, as well as a quick reference guide that you can refer to on the go. These will guide you through the easy steps to maintaining correct alignment with Pro-Flex LP Align. We also recommend that you practice with your prosthetic technician/professional to ensure your alignment is correct.

Want to know how to stay aligned? Check out our video or read the Quick Reference Guide.

WWW.OSSUR.COM

Heel height can be adjusted up to 7cm at the simple touch of a button

Full effective toe length for better stance and the use of sandals



Flexible toe accommodates adjustments for different heel heights

3-blade design provides dynamic energy return

