

Formfit® Post-Op Knee

Formfit Post-Op Knee provides range of motion control and immobilization following injury or surgery. Length adjustment is tool free and easy. The low profile ROM hinge can be adjusted with a simple pull, rotate and release motion. Wide range of flexion and extension settings allow you to choose the optimal range of motion.

Features and benefits

- Breakaway struts for easy, tool free length adjustment
- OneTouch™ easy to adjust ROM hinge with lock-out capability from -10° to 30° (10° incremental)
- Extension control from -10° to 90° (10° incremental)
- Flexion control from 0° to 120° (10° incremental)
- Available in 3 models, 3 lengths, and cool or full foam

Indications for use

- Post-op/post-injury ROM control or immobilization. Indicated for ligament and fracture-related surgeries and when controlled range of motion is necessary

Required Measurements

- Thigh circumference 15 cm (6") above mid patella

Ordering Information

Description	Thigh Circumference	Cool	Full Foam
Knee PO Cool Reg 24in DL	UP TO 27" (69cm)	FPOK403100	FPOK403200
Knee PO CoolXL 24in DL	27" – 35" (69-89cm)	FPOK403102	FPOK403202
Knee PO Cool Reg 28in DL	UP TO 27" (69cm)	FPOK403101	FPOK403201
Knee PO Cool XL 28in DL	27" – 35" (69-89cm)	FPOK403103	FPOK403203
Knee PO Full Reg 16in DL	UP TO 27" (69cm)	-	FPOK403300
Knee PO Full XL 16in DL	27" – 35" (69-89cm)	-	FPOK403301

Soft Goods

Description	Cool	Full Foam
Replacement Softgood, Regular	403175	403275
Replacement Softgood, XL	403176	403277
Replacement Softgood, Regular Long	403276	403278

Accessories

Description	
ROM stop locks (set of 20 lock clips)	B-8380005
Strap Reduction Kit Cool	403199
Strap Reduction Kit Full Foam	403299

