

MESH TOP POST-OP SHOE

INDICATIONS FOR USE

- Post-operative use or minor trauma wounds
- Geriatric or insensitive foot

FEATURES AND BENEFITS

- Breathable mesh upper with reinforced heel counter and padded heel collar provides increased stability overall to the rear of the foot while reducing slippage
- Semi-rigid sole provides superior protection and control after osteotomies
- Extended strap and dual side closure system accommodates extra bulky dressings



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Shoe Size UK		Shoe Size EU			
	Male	Female	Male	Female	Male	Female
Small	3.5 - 5.5	2.5 - 3.5	35.5 - 38	35 - 36	17013	17003
Medium	6 - 8	4 - 5.5	38.5 - 42	36.5 - 38.5	17015	17005
Large	8.5 - 10.5	6 +	42.5 - 45	39 +	17017	17007
X-Large	11 - 13.5		45.5 - 48.5		17018	