



## LIVINGSKIN™

A passive functional prosthesis, hand crafted from high definition silicone to match exact skin tone and appearance.

Össur offers passive Livingskin™ prostheses for:

partial fingers



whole fingers



partial hands



hands



below elbow and above elbow





LivingSkin™ is appropriate for anyone with any level of upper limb absence who would prefer a life-like and functionally supportive prosthesis.

Manufactured from versatile, durable and virtually stain-proof silicone.

Adjustable armatures may be fitted to allow repositioning of each finger for specific tasks.

Realistic in appearance - matching of exact skin tone, addition of freckles, hair, scars and tattoos.

Ability to assist in Activities of Daily Living such as pushing, pulling, stabilizing, supporting, light grasping and typing.

FOLLOW ÖSSUR ON



USA (800) 233-6263  
[WWW.OSSUR.COM](http://WWW.OSSUR.COM)

CANADA (800) 663-5982  
[WWW.OSSUR.CA](http://WWW.OSSUR.CA)

